

TECHNICAL BULLETIN

⚠ **Caution:** It is important to completely prime the pump tubing. Failure to do so may prevent the pump from infusing.

Use Aseptic Technique

The following steps should be followed when priming the Halyard elastomeric pumps.

- Open clamp.
- Remove tubing cap to start priming (up to 15 minutes).
- When all air has been removed from the entire tubing and fluid flow is observed at end of the distal luer, the administration set is primed.
- Close clamp and replace tubing cap until ready for use.

If the administration set does not prime, follow these steps:

1. Attach a luer adapter or stopcock to the distal luer.
2. Attach a small syringe (10 ml preferred) to the other side of the adapter and pull back on the syringe to create suction.
3. Continue to create suction until all air is removed from the tubing and fluid flow is observed from the distal luer. Repeat as necessary.
4. Disconnect syringe and luer adaptor or stopcock, and observe pump for complete priming.
5. If this does not work, check to see if something else is impeding flow, such as medication precipitate, clamp is closed or kinked tubing.

NOTE: The above procedures are not applicable to ON-Q* pumps with Select-A-Flow* and/or ONDEMAND* devices.

There are inherent risks in all medical devices. Please refer to the product labeling for **Indications, Cautions, Warnings and Contraindications**. Failure to follow the product labeling could directly impact patient safety. Physician is responsible for prescribing and administering medications per instructions provided by the drug manufacturer. Refer to www.halyardhealth.com for additional product safety **Technical Bulletins**.

Please contact the Clinical Services Department at **800-444-2728** or **949-923-2400**
if you have any questions regarding this information.
