

## Non-Adherent Dry Dressings

# Adhesive Island Dressings

- General purpose bordered gauze dressing
- Easy to use all-in-one dressing is conformable and a gentle alternative to gauze and tape
- Non-adherent contact layer with absorptive pad
- Soft, absorbent pad allows light to moderate fluid absorption without adhering to the wound or disrupting healing tissue
- Non-woven, breathable and non-irritating adhesive border with rounded edges secures the dressing in place and protects the wound
- Split release liner for easy, one-step application
- Ideal for incision sites and other wounds as a primary or secondary dressing
- Available in a range of sizes
- Available in standard as well as plus option for increased fluid absorption

**NOT MADE WITH**  
natural rubber latex



Single Use Only

hypoallergenic

Sterile		
Catalog #	Description	Quantity
IDS61001	2" x 3" (5.1 x 7.6 cm)	50 EA/BX, 8 BX/CA
IDS61002	4" x 4" (10.2 x 10.2 cm)	5 EA/BX, 8 BX/CA
IDS61003	4" x 5" (10.2 x 12.7 cm)	25 EA/BX, 8 BX/CA
IDS61004	4" x 8" (10.2 x 20.3 cm)	25 EA/BX, 4 BX/CA
IDS61005	4" x 10" (10.2 x 25.4 cm)	25 EA/BX, 4 BX/CA
IDS61006	4" x 14" (10.2 x 35.6 cm)	25 EA/BX, 2 BX/CA
IDS61007	6" x 6" (15.2 x 15.2 cm)	25 EA/BX, 4 BX/CA
IDS62001	4" x 6" (10.2 x 15.2 cm)	25 EA/BX, 4 BX/CA
IDS62002	6" x 7" (15.2 x 17.8 cm)	25 EA/BX, 4 BX/CA

