


# INTRODUCING HALYARD\* BORDERED FOAM DRESSINGS

**PERFORMANCE. COMFORT. SAVINGS. TREAT WITH CONFIDENCE.**

**HALYARD\* multi-layered bordered silicone foam dressing** made with a superabsorbent foam layer, protective film backing and gentle silicone adhesive.

- For management of moderate-to-heavy drainage, partial and full thickness wounds
- Multi-layer design to protect and absorb exudate from the wound
- Promote a moist wound-healing environment
- May remain in place for up to 7 days
- May be used as part of a pressure injury prevention protocol



**OUR NEW WOUND CARE LINE IS PART OF THE HALYARD\* PORTFOLIO**



INFECTION CONTROL



INCONTINENCE




SURGICAL SOLUTIONS



ADVANCED WOUND CARE

## Bordered Foam Dressings

Scientific testing showed HALYARD\* Bordered Silicone Foam Dressing performs as well as the market leading foam dressings in both conformability and absorption.



Source: HALYARD\* data on file

PRODUCT CODE	DESCRIPTION	CASE QTY.
50084	HALYARD* Multi-Layer Bordered Silicone Foam Dressing, 3 in. x 3 in. (7.6cm x 7.6cm)	100
50085	HALYARD* Multi-Layer Bordered Silicone Foam Dressing, 4 in. x 4 in. (10cm x 10cm)	100
55571	HALYARD* Multi-Layer Bordered Silicone Foam Dressing, 5.75 in. x 5.75 in. (14 cm x 14 cm)	100
50086	HALYARD* Multi-Layer Bordered Silicone Foam Dressing, 6 in. x 6 in. (15cm x 15cm)	100
55572	HALYARD* Multi-Layer Bordered Silicone Foam Dressing, 6 in. x 8 in. (15 cm x 20 cm)	100
50087	HALYARD* Multi-Layer Bordered Silicone Foam Dressing, Sacral, 7 in. x 6.5 in. (17.7cm x 16.5cm)	40
50088	HALYARD* Multi-Layer Bordered Silicone Foam Dressing, Sacral, 9 in. x 8 in. (22.8cm x 20cm)	25
55570	HALYARD* Multi-Layer Bordered Silicone Foam Dressing, Heel, 9 in. x 9 in. (22 cm x 22 cm)	25

Partner with O&M Halyard to stay in the forefront of patient care.  
For more information, please visit: [www.halyardhealth.com](http://www.halyardhealth.com) or contact your  
O&M Halyard representative **1-844-HALYARD** (1-844-425-9273).