

ASSURED FACIAL PROTECTION CAN BE AS EASY AS 1,2,3

From standard and fluid resistant masks, to particulate filter respirators, eyewear, and face shields, HALYARD* offers a broad range of styles, features, and levels of protection. HALYARD* face masks are the choice of more healthcare professionals where comfort and protection are needed in a variety of tasks and procedures.

HALYARD* LEADS THE WAY IN ASTM-RATED PROTECTION

HALYARD*, a world leader in face masks, offers products that meet the ASTM F2100 standard at all levels. This helps you make sure you're protected appropriately for the fluid risk that may be present. You can, at a glance, pick the right mask for the task. It's as easy as 1-2-3.

ASTM International is a voluntary standards development organization whose work covers 15 industries, ranging from energy to chemicals to medical devices. ASTM F2100 is recognized by the FDA as the standard for medical face masks.



ASTM F2100: THE STANDARD FOR FLUID-RESISTANT PROTECTION

ASTM has issued guidelines for performance of medical face masks, including a requirement for fluid resistance and particulate filtration for all masks. ASTM F2100 standards also call for mask packaging to be clearly labeled with the level of protection at one of three levels: 1, 2, or 3. This standard reinforces the importance of labelling to identify protection level, based on fluid resistance, bacterial filtration, particulate filtration, and breathability. ASTM F2100 is recognized by the FDA as the consensus standard for medical face masks.

		LEVEL 1	LEVEL 2	LEVEL 3
Test Method	Test Description	ASTM F2100-19	ASTM F2100-19	ASTM F2100-19
ASTM F1862	Fluid Resistance with synthetic blood, in mm Hg	80 mm Hg	120 mm Hg	160 mm Hg
EN14683 Delta P	Delta P Differential pressure, mm H ₂ O/cm ²	<5.0 mm H ₂ O	<6.0 mm H ₂ O	<6.0 mm H ₂ O
ASTM F2101	Bacterial Filtration Efficiency (BFE), %	≥ 95%	≥ 98%	≥ 98%
ASTM F2299	Submicron Particle Filtration (PFE) @ 0.1 μ, %	≥ 95%	≥ 98%	≥ 98%
16 CFR Part 1610	Flammability	Class 1	Class 1	Class 1

THE RIGHT MASK FOR THE RIGHT TASK



Procedure/Earloop Mask

Masks with earloops designed for easy donning typically used in patient care areas



Surgical Ties

Masks with ties designed for a more custom, secure fit typically used in the OR and other "red lined" environments



Respirators

Respiratory protection designed to guard against airborne hazards

LEVEL 1	LEVEL 2	LEVEL 3
 For general purpose medical procedures where the wearer is not at risk of blood or bodily fluid splash or to protect staff and/or the patient from droplet exposure to microorganisms (e.g. patient with upper respiratory tract infection visits GP).	 For emergency departments, dentistry, changing dressings on small or healing wounds where minimal blood droplet exposure may possibly occur (e.g. endoscopy procedures).	 For all surgical procedures, major trauma first aid or in any area where the health care worker is at risk of blood or bodily fluid splash (e.g. orthopaedic, cardio-vascular procedures).

PROCEDURE MASKS WITH EARLOOPS

STANDARD PROCEDURE MASKS

Item #	Product Description	Items/Box	Boxes/Case	Case Count
00152	Standard Cone Mask, With Headband, Not Earloops, Blue	50	6	300
47080	Procedure Mask, SO SOFT* Earloops, Blue	50	10	500
47117	Procedure Mask, SO SOFT* Earloops, Yellow	50	10	500
47295	TEDDY BEAR* Procedure Mask with SO SOFT* Earloops, Teddy Bear Print	50	10	500
62356	THE LITE ONE* Procedure Mask, SO SOFT* Earloops, Blue	50	10	500
62363	Anti-Fog Procedure Mask with SO SOFT* Lining, White	50	10	500

LEVEL 1 PROCEDURE MASKS

Item #	Product Description	Items/Box	Boxes/Case	Case Count
25867	FLUIDSHIELD* Procedure Mask, SO SOFT* Lining and Earloops, Yellow	50	10	500
25868	FLUIDSHIELD* Procedure Mask, SO SOFT* Lining and Earloops, Lavender	50	10	500
25869	FLUIDSHIELD* Procedure Mask, SO SOFT* Lining and Earloops, Blue	50	10	500
41802	FLUIDSHIELD* Procedure Mask, Fog-Free Foam Strip, SO SOFT* Lining and Earloops, White	50	10	500
47566	Procedure Mask, White	50	10	500
47567	Procedure Mask, Yellow	50	10	500
47568	Procedure Mask, Blue	50	10	500
48388	Procedure Mask, SO SOFT* Earloops, Yellow	50	10	500
49513	Level 1 HALYARD* Fog-Free Earloop Mask, Black	50	10	500
57578	Procedure Mask, Pink	50	10	500

LEVEL 2 FLUIDSHIELD* PROCEDURE MASKS

Item #	Product Description	Items/Box	Boxes/Case	Case Count
47573	Procedure Mask, Blue	50	10	500
62115	Level 2 FLUIDSHIELD* Procedure Mask, SO SOFT* Earloops, Blue	50	10	500
62116	Level 2 FLUIDSHIELD* Procedure Mask, WrapAround Visor, Fog-Free Foam Strip, SO SOFT* Earloops, Blue	25	4	100

LEVEL 3 FLUIDSHIELD* PROCEDURE MASKS

Item #	Product Description	Items/Box	Boxes/Case	Case Count
00146	Level 3 FLUIDSHIELD* Procedure Mask, WrapAround Visor, Fog-Free Foam Strip, SO SOFT* Lining and Earloops, Blue	25	4	100
00148	Level 3 FLUIDSHIELD* Procedure Mask, Fog-Free Foam Strip, SO SOFT* Lining and Earloops, Blue	40	10	400
28797	Level 3 FLUIDSHIELD* Procedure Mask, Fog-Free Foam Strip, SO SOFT* Earloops, Orange	40	10	400
28800	Level 3 FLUIDSHIELD* Procedure Mask, WrapAround Visor, Fog-Free Foam Strip, SO SOFT* Earloops, Orange	25	4	100
47107	Level 3 FLUIDSHIELD* Fog-Free Procedure Mask, SO SOFT* Lining & Earloops, Orange	40	10	400
47137	Level 3 FLUIDSHIELD* Procedure Mask, Visor, Fog-Free Foam Strip, SO SOFT* Lining and Earloops, Orange	25	4	100
47147	Level 3 FLUIDSHIELD* Fog-Free Procedure Mask, Anti-Glare WrapAround Visor, SO SOFT* Lining and Earloops, Orange	25	4	100
47298	Level 3 FLUIDSHIELD* Procedure Mask, WrapAround Visor, Fog-Free Foam Strip, SO SOFT* Lining and Earloops, Teddy Bear	25	4	100
47569	Procedure Mask, Anti-Fog Foam Strip, Aqua	50	10	500
49514	Level 3 HALYARD* Fog-Free Earloop Mask, Black	50	10	500

PEDIATRIC PATIENT MASK

Item #	Product Description	Items/Box	Boxes/Case	Case Count
32856	Procedure Mask, SO SOFT* Lining, Child, Disney® Characters, Ages 4-12	75	10	750

SURGICAL MASKS WITH TIES

STANDARD SURGICAL MASKS

Item #	Product Description	Items/Box	Boxes/Case	Case Count
48100	THE LITE ONE* Surgical Mask, Blue	50	6	300
48390	Surgical Mask with SO SOFT* Lining, White	50	6	300
49214	Anti-Fog Surgical Mask, Blue	50	6	300

STANDARD SPECIALTY SURGICAL MASKS

Item #	Product Description	Items/Box	Boxes/Case	Case Count
37525	DUCKBILL* Surgical Mask, Green	50	6	300
47500	SOFT TOUCH* II Surgical Mask, Horizontal Ties, Blue	50	6	300
47700	Pouch Style Mask with SO SOFT* Lining, Headband, Not Ties, Blue	50	6	300
48220	DUCKBILL* Surgical Mask, Blue	50	6	300
49215	Anti-Fog Surgical Mask, with DERMA-TOUCH* Tape, Green	50	6	300
49216	DUCKBILL* Anti-Fog Surgical Mask, Blue	50	6	300
49235	Anti-Fog Surgical Mask, with Film Strip, Green	50	6	300
49310	ULTI-MASK* Anti-Fog Surgical Mask, White	50	6	300

LEVEL 1 FLUIDSHIELD* SURGICAL MASKS

Item #	Product Description	Items/Box	Boxes/Case	Case Count
28802	Surgical Mask with SO SOFT* Lining, Blue	50	6	300
28806	Anti-Fog Surgical Mask with SO SOFT* Lining, Green	50	6	300
41805	Surgical Mask with SO SOFT* Lining, White	50	6	300
47571	Level 1 Surgical Mask, Softbill, Blue	50	6	300
47578	Surgical Mask, White	50	6	300

LEVEL 1 FLUIDSHIELD* SPECIALTY SURGICAL MASKS

Item #	Product Description	Items/Box	Boxes/Case	Case Count
39118	Surgical Mask, with SO SOFT* Lining, High Clarity WrapAround Visor, Horizontal Ties, Blue	25	4	100
47570	Surgical Mask, Horizontal Ties, Blue	50	6	300
47572	Surgical Mask Blue	50	6	300
49325	Level 1 Anti-Fog Surgical Mask, Derma-Touch Tape, Green	50	6	300

LEVEL 2 FLUIDSHIELD* SURGICAL MASKS

Item #	Product Description	Items/Box	Boxes/Case	Case Count
62113	Anti-Fog Surgical Mask, Blue with Orange Print	50	6	300
62114	Anti-Fog Surgical Mask, WrapAround Visor, Blue with Orange Print	25	4	100

LEVEL 2 FLUIDSHIELD* SPECIALTY SURGICAL MASKS

Item #	Product Description	Items/Box	Boxes/Case	Case Count
39123	Expanded Chamber Surgical Mask, with SO SOFT* Lining, Blue/White	50	6	300
39124	Expanded Chamber Surgical Mask, with SO SOFT* Lining, High Clarity WrapAround Visor, Blue/White	25	4	100

LEVEL 3 FLUIDSHIELD* SURGICAL MASKS

Item #	Product Description	Items/Box	Boxes/Case	Case Count
28804	Anti-Fog Surgical Mask, Anti-Glare WrapAround Visor, Orange	25	4	100
47574	Surgical Mask, Aqua	50	6	300
47575	Anti-Fog Surgical Mask, Aqua	50	6	300
48207	Anti-Fog Surgical Mask with SO SOFT* Lining, Orange	50	6	300
48237	Anti-Fog Surgical Mask with SO SOFT* Lining, Visor, Orange	25	4	100
48247	Anti-Fog Surgical Mask with SO SOFT* Lining, Anti-Glare WrapAround Visor, Orange	25	4	100
48297	TEDDY BEAR* Surgical Mask with SO SOFT* Lining, Teddy Bear Print	50	6	300
49326	Level 3 Anti-Fog Surgical Mask, Derma-Touch Tape, Green	50	6	300

EYEWEAR AND FACE SHIELDS

SAFEVIEW* EYEWEAR

Item #	Description	Items/Box	Items/Case
SV50A	Eyewear Assembled Glasses, Assorted Colors	5	50
50005	SAFEVIEW CURVE* Comfort Eyewear Dispensing Tower	4	100
48418	SAFEVIEW CURVE* Wallbracket	1	1

SAFEVIEW* EYEWEAR COMPONENTS

Item #	Description	Items/Bag	Items/Case
SV100F	Eyewear Frames, Assorted Colors	10	100
SV250L	Eyewear Replacement Lenses	10	250



HALYARD* FACE SHIELDS

Item #	Description	Items/Bags	Bags/Case
41204	Face Shield, Anti-Fog Film with Foam Band, Headband, 9.5 inch Long	40	1
41205	Face Shield, Anti-Fog Film with Foam Band, Headband, 7.5 inch Long	40	1

N95 PARTICULATE FILTER RESPIRATORS

FLUIDSHIELD* N95 PARTICULATE FILTER RESPIRATORS

Item #	Product Description	Items/Box	Boxes/Case	Case Count
46727	N95 Particulate Filter Respirator and Surgical Mask with SO SOFT* Lining, Regular, Orange	35	6	210
46827	N95 Particulate Filter Respirator and Surgical Mask with SO SOFT* Lining, Small, Orange	35	6	210
46767	Particulate Filter Respirator and Surgical Mask, Safety Seal Film , SO SOFT* Lining, Regular, Orange	35	6	210
46867	Particulate Filter Respirator and Surgical Mask, Safety Seal Film , SO SOFT* Lining, Small, Orange	35	6	210
49299	HALYARD* N95 Surgical Cup Respirator Masks, Small (smaller to narrower faces), Blue	15	12	180
49300	HALYARD* N95 Surgical Cup Respirator Masks, Regular (average to large faces), Blue	15	12	180
62126	N95 Particulate Filter Respirator and Surgical Mask with SO SOFT* Lining, Regular, White	50	6	300

N95 PARTICULATE FILTER RESPIRATORS ACCESSORIES

Item #	Product Description	Items/Box	Boxes/Case	Case Count
47950	Qualitative Fit Test Kit	1	1	1
99003	Saccharin Solutions for Qualitative Fit Testing	2	1	2
99009	Bitrex Solutions for Qualitative Fit Testing	2	1	2
99010	Hood Assemblies for Qualitative Fit Testing, Yellow	10	1	10
99015	Hood Assemblies with Bouffant Cap for Qualitative Fit Testing, Clear	2	1	2
99018	Nebulizers for Qualitative Fit Testing	2	1	2

To place orders online, access the customer portal (www.hyhcustomerportal.com) and select the Order Online tab. To speak to a an HALYARD* Customer Service Representative, contact 844-HALYARD (844-425-9273) or email MyHalyardOrders@owens-minor.com.

Please visit our product catalog for our full listing:
<https://products.halyardhealth.com/infection-prevention/facial-respiratory-protection.html>

For more information, please visit: www.halyardhealth.com

Call 1-844-HALYARD (1-844-425-9273) in the United States and Canada.

www.halyardhealth.com



HALYARD* is an Owens & Minor brand

Disney Elements ©DISNEY

*Registered Trademark or Trademark of O&M Halyard or its affiliates. ©2018, 2024, 2025. All rights reserved. COPY-05089 008906