



Magic in Care*



Where Comfort Meets Protection

HALYARD* LAVENDER* Nitrile Exam Gloves for Disney makes protection fun! These exam gloves are thinner, lighter, and economical, yet they retain the protective properties required for excellent performance across a wide range of healthcare tasks. So, they provide the nitrile glove protection you can trust.



Unique dispenser color and Disney character for each glove size!

Features:

- Colorful, child-friendly dispenser design inspired by the Disney Frozen characters
- Powder-free
- Textured fingertips
- Beaded cuffs
- Low Dermatitis Potential (LDP) tested
- Tested with 42 chemo & hazardous drugs
- Tested with fentanyl & gastric acid
- Not made with natural rubber latex

Item	Size	Packaging
50067	X-Small	250 eaches/box 10 boxes/case
50068	Small	250 eaches/box 10 boxes/case
50069	Medium	250 eaches/box 10 boxes/case
50070	Large	250 eaches/box 10 boxes/case
50071	X-Large	230 eaches/box 10 boxes/case





Magic in Care*



Product Specifications

Gauge Thickness Measurements	MM	MIL
Middle Finger:	.08	3.1
Palm:	.06	2.3
Cuff:	.05	2.0
Average Length	242mm	(9.5")

Physical Properties

Before Aging

Tensile Strength	23 MPa
Ultimate Elongation	600%

After Aging

Tensile Strength	28 MPa
Ultimate Elongation	550%

All physical properties and product specs are target measurements; target MIL measurements may differ due to rounding.

Quality Standards

Exceeds current ASTM Standards for critical defects, tensile strength, and puncture resistance.

No evidence of dermal irritation or contact sensitization as measured by Primary Skin Irritation, Repeat Challenge Sensitization.

Other Benefits

- LAVENDER* Nitrile Exam Gloves have an AQL of 1.0 for critical defects. This means that LAVENDER* Exam Gloves are significantly less likely to have critical defects than exam gloves that just meet the ASTM 2.5 AQL requirement.

Ordering Information

To place orders online, access the customer portal (www.hyhcustomerportal.com) and select the Order Online tab. To speak to a HALYARD* Customer Service Representative, contact 844-HALYARD (844-425-9273) or email MyHalyardOrders@owens-minor.com