

FACIAL PROTECTION



Procedure Masks with Earloops

COMFORT MEETS PROTECTION

Select the Correct Mask for the Task

Procedure masks with earloops are designed for easy donning and are typically used in patient care areas.

HALYARD*, a world leader in face masks products that meet the ASTM F2100 standard at all levels. This helps make sure you're protected appropriately for the fluid risk that may be present. You can, at a glance, pick the right mask for the task. It's as easy as 1-2-3. ASTM International is a voluntary standards development organization whose work covers 15 industries, ranging from energy to chemicals to medical devices. ASTM F2100 is recognized by the FDA as the standard for all medical face masks.

ASTM has issued guidelines for performance of medical face masks, including a requirement for fluid resistance and particulate filtration for all masks. ASTM F2100 standards also call for mask packaging to be clearly labeled with the level of protection at one of three levels: 1, 2, or 3. This standard reinforces the importance of labelling to identify protection level, based on fluid resistance, bacterial filtration, particulate filtration, and breathability.

Features

- All levels of protection + non-rated
- SO SOFT* lining and earloops
- Adjustable nosepieces
- Fog-free features
- WrapAround and Anti-glare Visors
- Different shapes and styles
- Variety of colors/designs



FACIAL PROTECTION



Procedure Masks with Earloops

STANDARD PROCEDURE MASKS

Item #	Product Description	Items/Box	Boxes/Case	Case Count
00152	Standard Cone Mask, With Headband, Not Earloops, Blue	50	6	300
47080	Procedure Mask, SO SOFT* Earloops, Blue	50	10	500
47117	Procedure Mask, SO SOFT* Earloops, Yellow	50	10	500
47295	TEDDY BEAR* Procedure Mask with SO SOFT* Earloops, Teddy Bear Print	50	10	500
62356	THE LITE ONE* Procedure Mask, SO SOFT* Earloops, Blue	50	10	500
62363	Anti-Fog Procedure Mask with SO SOFT* Lining, White	50	10	500

LEVEL 1 PROCEDURE MASKS

Item #	Product Description	Items/Box	Boxes/Case	Case Count
25867	FLUIDSHIELD* Procedure Mask, SO SOFT* Lining and Earloops, Yellow	50	10	500
25868	FLUIDSHIELD* Procedure Mask, SO SOFT* Lining and Earloops, Lavender	50	10	500
25869	FLUIDSHIELD* Procedure Mask, SO SOFT* Lining and Earloops, Blue	50	10	500
41802	FLUIDSHIELD* Procedure Mask, Fog-Free Foam Strip, SO SOFT* Lining and Earloops, White	50	10	500
47566	Procedure Mask, White	50	10	500
47567	Procedure Mask, Yellow	50	10	500
47568	Procedure Mask, Blue	50	10	500
48388	Procedure Mask, SO SOFT* Earloops, Yellow	50	10	500
49513	Level 1 HALYARD* Fog-Free Earloop Mask, Black	50	10	500
57578	Procedure Mask, Pink	50	10	500

LEVEL 2 FLUIDSHIELD* PROCEDURE MASKS

Item #	Product Description	Items/Box	Boxes/Case	Case Count
47573	Procedure Mask, Blue	50	10	500
62115	Level 2 FLUIDSHIELD* Procedure Mask, SO SOFT* Earloops, Blue	50	10	500
62116	Level 2 FLUIDSHIELD* Procedure Mask, WrapAround Visor, Fog-Free Foam Strip, SO SOFT* Earloops, Blue	25	4	100

LEVEL 3 FLUIDSHIELD* PROCEDURE MASKS

Item #	Product Description	Items/Box	Boxes/Case	Case Count
00146	Level 3 FLUIDSHIELD* Procedure Mask, WrapAround Visor, Fog-Free Foam Strip, SO SOFT* Lining and Earloops, Blue	25	4	100
00148	Level 3 FLUIDSHIELD* Procedure Mask, Fog-Free Foam Strip, SO SOFT* Lining and Earloops, Blue	40	10	400
28797	Level 3 FLUIDSHIELD* Procedure Mask, Fog-Free Foam Strip, SO SOFT* Earloops, Orange	40	10	400
28800	Level 3 FLUIDSHIELD* Procedure Mask, WrapAround Visor, Fog-Free Foam Strip, SO SOFT* Earloops, Orange	25	4	100
47107	Level 3 FLUIDSHIELD* Fog-Free Procedure Mask, SO SOFT* Lining & Earloops, Orange	40	10	400
47137	Level 3 FLUIDSHIELD* Procedure Mask, Visor, Fog-Free Foam Strip, SO SOFT* Lining and Earloops, Orange	25	4	100
47147	Level 3 FLUIDSHIELD* Fog-Free Procedure Mask, Anti-Glare WrapAround Visor, SO SOFT* Lining and Earloops, Orange	25	4	100
47298	Level 3 FLUIDSHIELD* Procedure Mask, WrapAround Visor, Fog-Free Foam Strip, SO SOFT* Lining and Earloops, Teddy Bear	25	4	100
47569	Procedure Mask, Anti-Fog Foam Strip, Aqua	50	10	500
49514	Level 3 HALYARD* Fog-Free Earloop Mask, Black	50	10	500

PEDIATRIC PATIENT MASK

Item #	Product Description	Items/Box	Boxes/Case	Case Count
32856	Procedure Mask, SO SOFT* Lining, Child, Disney® Characters, Ages 4-12	75	10	750

FACIAL PROTECTION




Procedure Masks with Earloops

ASTM F2100: THE STANDARD FOR FLUID-RESISTANT PROTECTION

ASTM has issued guidelines for performance of medical face masks, including a requirement for fluid resistance and particulate filtration for all masks. ASTM F2100 standards also call for mask packaging to be clearly labeled with the level of protection at one of three levels: 1, 2, or 3. This standard reinforces the importance of labelling to identify protection level, based on fluid resistance, bacterial filtration, particulate filtration, and breathability. ASTM F2100 is recognized by the FDA as the consensus standard for medical face masks.

Test Method	Test Description	LEVEL 1	LEVEL 2	LEVEL 3
		ASTM F2100-19	ASTM F2100-19	ASTM F2100-19
ASTM F1862	Fluid Resistance with synthetic blood, in mm Hg	80 mm Hg	120 mm Hg	160 mm Hg
EN14683 Delta P	Delta P Differential pressure, mm H ₂ O/cm ²	<5.0 mm H ₂ O	<6.0 mm H ₂ O	<6.0 mm H ₂ O
ASTM F2101	Bacterial Filtration Efficiency (BFE), %	≥ 95%	≥ 98%	≥ 98%
ASTM F2299	Submicron Particle Filtration (PFE) @ 0.1 µ, %	≥ 95%	≥ 98%	≥ 98%
16 CFR Part 1610	Flammability	Class 1	Class 1	Class 1

THE RIGHT MASK FOR THE RIGHT TASK



Procedure/Earloop Mask

Masks with earloops designed for easy donning typically used in patient care areas

Surgical Ties

Masks with ties designed for a more custom, secure fit typically used in the OR and other "red lined" environments

Respirators

Respiratory protection designed to guard against airborne hazards

Ordering Information

To place orders online, access the customer portal (www.hyhcustomerportal.com) and select the Order Online tab. To speak to a HALYARD* Customer Service Representative, contact **844-HALYARD (844-425-9273)** or email MyHalyardOrders@owens-minor.com