



Magic in Care*



Where Comfort Meets Protection

HALYARD* **STERLING SG* Nitrile Exam Gloves** for Disney makes protection fun! This exam glove offers a comfortable fit and feel and is the worry-free choice – designed for use by virtually everyone, everywhere. If you want quality performance in a nitrile glove and want to standardize on one exam glove for virtually all your needs, facility-wide, **STERLING SG* Nitrile** is an obvious choice.



Unique dispenser color
and Disney character
for each glove size!

Features:

- Colorful, child-friendly dispenser design inspired by the Disney characters
- Powder-free
- Textured fingertips
- Beaded cuffs
- Low Dermatitis Potential (LDP) tested
- Tested with 52 chemo & hazardous drugs
- Tested with fentanyl & gastric acid
- Not made with natural rubber latex

Item	Size	Packaging
50238	X-Small	250 eaches/box 10 boxes/case
50239	Small	250 eaches/box 10 boxes/case
50240	Medium	250 eaches/box 10 boxes/case
50241	Large	250 eaches/box 10 boxes/case
50242	X-Large	230 eaches/box 10 boxes/case





Magic in Care*



Product Specifications

Gauge Thickness Measurements	MM	MIL
Middle Finger:	.09	3.5
Palm:	.08	3.1
Cuff:	.06	2.4
Average Length	242mm	(9.5")

Physical Properties

Before Aging

Tensile Strength	42 MPa
Ultimate Elongation	650%

After Aging

Tensile Strength	38 MPa
Ultimate Elongation	550%

All physical properties and product specs are target measurements; target MIL measurements may differ due to rounding.

Quality Standards

Exceeds current ASTM D 6319: Standard Specification for Nitrile Examination Gloves for Medical Application.

No evidence of dermal irritation or contact sensitization as measured by Primary Skin Irritation, Repeat Challenge Sensitization.

Note: All specifications are subject to change without notice.

Other Benefits

- SMARTPULL* Box Opening – significantly reducing the incidence of getting more than one glove while dispensing.
- Offers a comfortable fit and feel.
- STERLING SG* Nitrile Exam Gloves have an AQL of 1.0 for critical defects. This means that STERLING SG* Exam Gloves are significantly less likely to have critical defects than exam gloves that just meet the ASTM 2.5 AQL requirement.

Ordering Information

To place orders online, access the customer portal (www.hyhcustomerportal.com) and select the Order Online tab. To speak to a HALYARD* Customer Service Representative, contact 844-HALYARD (844-425-9273) or email MyHalyardOrders@owens-minor.com