

SELECT THE CORRECT MASK FOR THE TASK

When it comes to selecting facial protection in a healthcare setting it is important to understand where the mask or respirator will be worn and potential exposure to blood, bodily fluids and surgical smoke. As compared to industrial N95s, surgical N95 respirators are constructed and tested to meet requirements for breathability, fluid resistance, flammability and particle filtration.

Not All N95 Respirators are Created Equal

- **Increased Comfort**
 - Focused on Breathability and Fit
- **Highest Level of Fluid Resistance**
- **Reliable Supply**
 - Benefits of a Vertically Integrated Supply Chain



FLUIDSHIELD* N95 PARTICULATE FILTER RESPIRATORS

Item #	Product Description	Items/Box	Boxes/Case	Case Count
46727	HALYARD* FLUIDSHIELD* 3 N95 Particulate Filter Respirator and Surgical Mask with SO SOFT* Lining, Regular, Orange	35	6	210
46827	HALYARD* FLUIDSHIELD* 3 N95 Particulate Filter Respirator and Surgical Mask with SO SOFT* Lining, Small, Orange	35	6	210
46767	HALYARD* FLUIDSHIELD* 3 Particulate Filter Respirator and Surgical Mask, Safety Seal Film , SO SOFT* Lining, Regular, Orange	35	6	210
46867	HALYARD* FLUIDSHIELD* 3 Particulate Filter Respirator and Surgical Mask, Safety Seal Film , SO SOFT* Lining, Small, Orange	35	6	210
49299	HALYARD* N95 Surgical Cup Respirator Masks, Small (smaller to narrower faces), Blue	15	12	180
49300	HALYARD* N95 Surgical Cup Respirator Masks, Regular (average to large faces), Blue	15	12	180
62126	HALYARD* FLUIDSHIELD* 2 N95 Particulate Filter Respirator and Surgical Mask with SO SOFT* Lining, Regular, White	50	6	300

TRUSTED IN HOSPITALS
FOR OVER TWENTY YEARS



PROTECTION THAT
EXCEEDS NIOSH
STANDARDS



REGULAR & SMALL SIZES
WITH SAME ASTM FLUID
RESISTANCE



VERTICALLY
INTEGRATED RAW
MATERIALS



EASY TO
SHAPE



LIGHTER
ON THE FACE



IDEAL FOR
STOCKPILING

Fit Test Accessories

Contains instructions and materials to conduct Fit Testing of N95 Respirators per requirements.

Features:

- 75–100 tests per kit
- Component Code
- Saccharin Threshold and Test Solutions 99003
- Yellow Tyvek Hood 99010
- Two nebulizers 99018 (one for threshold solution and one for test solution)
- Nebulizer Support Fitment and Tool (2)

Fit Test Accessories

Product Code	Description	Items/Case
	47950 Qualitative Fit Test Kit	1
	99003 Saccharin Solutions for Qualitative Fit Testing	2
	99009 BITREX® Solutions for Qualitative Fit Testing	2
	99010 Hood Assemblies for Qualitative Fit Testing, Yellow	10
	99015 Hood Assemblies with Bouffant Cap for Qualitative Fit Testing, Clear	2
	99018 Nebulizers for Qualitative Fit Testing	2

To place orders online, access the customer portal (www.hyhcustomerportal.com) and select the Order Online tab. To speak to a HALYARD* Customer Service Representative, contact **844-HALYARD** (844-425-9273) or email MyHalyardOrders@owens-minor.com