

FACIAL PROTECTION

Surgical Tie Masks



SELECT THE CORRECT MASK FOR THE TASK

Select the Correct Mask for the Task

Masks with ties are designed for a more custom, secure fit and are typically used in the OR and other "red-lined" environments. These masks provide a secure, customized fit that enhances comfort and ensures optimal protection.

HALYARD*, a world leader in face masks products that meet the ASTM F2100 standard at all levels. This helps make sure you're protected appropriately for the fluid risk that may be present. You can, at a glance, pick the right mask for the task. It's as easy as 1-2-3. ASTM International is a voluntary standards development organization whose work covers 15 industries, ranging from energy to chemicals to medical devices. ASTM F2100 is recognized by the FDA as the standard for all medical face masks.

ASTM has issued guidelines for performance of medical face masks, including a requirement for fluid resistance and particulate filtration for all masks. ASTM F2100 standards also call for mask packaging to be clearly labeled with the level of protection at one of three levels: 1, 2, or 3. This standard reinforces the importance of labelling to identify protection level, based on fluid resistance, bacterial filtration, particulate filtration, and breathability.

Features

- All levels of protection
- Pleated styles with SO SOFT* Inner lining
- Adjustable nosepieces
- Fog-free: Foam, Tape, Film
- WrapAround and Anti-glare Visors
- Expanded Chamber
- Variety of colors/designs
- Ties: Longer ties, Horizontal design



FACIAL PROTECTION



Surgical Tie Masks

STANDARD SURGICAL MASKS				
Item #	Product Description	Items/Box	Boxes/Case	Case Count
48100	THE LITE ONE* Surgical Mask, Blue	50	6	300
48390	Surgical Mask with SO SOFT* Lining, White	50	6	300
49214	Anti-Fog Surgical Mask, Blue	50	6	300
STANDARD SPECIALTY SURGICAL MASKS				
Item #	Product Description	Items/Box	Boxes/Case	Case Count
37525	DUCKBILL* Surgical Mask, Green	50	6	300
47500	SOFT TOUCH* II Surgical Mask, Horizontal Ties, Blue	50	6	300
47700	Pouch Style Mask with SO SOFT* Lining, Headband, Not Ties, Blue	50	6	300
48220	DUCKBILL* Surgical Mask, Blue	50	6	300
49215	Anti-Fog Surgical Mask, with DERMA-TOUCH* Tape, Green	50	6	300
49216	DUCKBILL* Anti-Fog Surgical Mask, Blue	50	6	300
49235	Anti-Fog Surgical Mask, with Film Strip, Green	50	6	300
49310	ULTI-MASK* Anti-Fog Surgical Mask, White	50	6	300
LEVEL 1 FLUIDSHIELD* SURGICAL MASKS				
Item #	Product Description	Items/Box	Boxes/Case	Case Count
28802	Surgical Mask with SO SOFT* Lining, Blue	50	6	300
28806	Anti-Fog Surgical Mask with SO SOFT* Lining, Green	50	6	300
41805	Surgical Mask with SO SOFT* Lining, White	50	6	300
47571	Level 1 Surgical Mask, Softbill, Blue	50	6	300
47578	Surgical Mask, White	50	6	300
LEVEL 1 FLUIDSHIELD* SPECIALTY SURGICAL MASKS				
Item #	Product Description	Items/Box	Boxes/Case	Case Count
39118	Surgical Mask, with SO SOFT* Lining, High Clarity WrapAround Visor, Horizontal Ties, Blue	25	4	100
47570	Surgical Mask, Horizontal Ties, Blue	50	6	300
47572	Surgical Mask Blue	50	6	300
49325	Level 1 Anti-Fog Surgical Mask, Derma-Touch Tape, Green NEW*	50	6	300
LEVEL 2 FLUIDSHIELD* SURGICAL MASKS				
Item #	Product Description	Items/Box	Boxes/Case	Case Count
62113	Anti-Fog Surgical Mask, Blue with Orange Print	50	6	300
62114	Anti-Fog Surgical Mask, WrapAround Visor, Blue with Orange Print	25	4	100
LEVEL 2 FLUIDSHIELD* SPECIALTY SURGICAL MASKS				
Item #	Product Description	Items/Box	Boxes/Case	Case Count
39123	Expanded Chamber Surgical Mask, with SO SOFT* Lining, Blue/White	50	6	300
39124	Expanded Chamber Surgical Mask, with SO SOFT* Lining, High Clarity WrapAround Visor, Blue/White	25	4	100
LEVEL 3 FLUIDSHIELD* SURGICAL MASKS				
Item #	Product Description	Items/Box	Boxes/Case	Case Count
28804	Anti-Fog Surgical Mask, Anti-Glare WrapAround Visor, Orange	25	4	100
47574	Surgical Mask, Aqua	50	6	300
47575	Anti-Fog Surgical Mask, Aqua	50	6	300
48207	Anti-Fog Surgical Mask with SO SOFT* Lining, Orange	50	6	300
48237	Anti-Fog Surgical Mask with SO SOFT* Lining, Visor, Orange	25	4	100
48247	Anti-Fog Surgical Mask with SO SOFT* Lining, Anti-Glare WrapAround Visor, Orange	25	4	100
48297	TEDDY BEAR* Surgical Mask with SO SOFT* Lining, Teddy Bear Print	50	6	300
49326	Level 3 Anti-Fog Surgical Mask, Derma-Touch Tape, Green NEW*	50	6	300

FACIAL PROTECTION

Surgical Tie Masks



ASTM F2100: THE STANDARD FOR FLUID-RESISTANT PROTECTION

ASTM has issued guidelines for performance of medical face masks, including a requirement for fluid resistance and particulate filtration for all masks. ASTM F2100 standards also call for mask packaging to be clearly labeled with the level of protection at one of three levels: 1, 2, or 3. This standard reinforces the importance of labelling to identify protection level, based on fluid resistance, bacterial filtration, particulate filtration, and breathability. ASTM F2100 is recognized by the FDA as the consensus standard for medical face masks.

Test Method	Test Description	LEVEL 1	LEVEL 2	LEVEL 3
		ASTM F2100-19	ASTM F2100-19	ASTM F2100-19
ASTM F1862	Fluid Resistance with synthetic blood, in mm Hg	80 mm Hg	120 mm Hg	160 mm Hg
EN14683 Delta P	Delta P Differential pressure, mm H ₂ O/cm ²	<5.0 mm H ₂ O/cm ²	<6.0 mm H ₂ O/cm ²	<6.0 mm H ₂ O/cm ²
ASTM F2101	Bacterial Filtration Efficiency (BFE), %	≥ 95%	≥ 98%	≥ 98%
ASTM F2299	Submicron Particle Filtration (PFE) @ 0.1 µ, %	≥ 95%	≥ 98%	≥ 98%
16 CFR Part 1610	Flammability	Class 1	Class 1	Class 1

THE RIGHT MASK FOR THE RIGHT TASK



Procedure/Earloop Mask

Masks with earloops designed for easy donning typically used in patient care areas



Surgical Ties

Masks with ties designed for a more custom, secure fit typically used in the OR and other "red lined" environments



Respirators

Respiratory protection designed to guard against airborne hazards

Ordering Information

To place orders online, access the customer portal (www.hyhcustomerportal.com) and select the Order Online tab. To speak to a HALYARD* Customer Service Representative, contact **844-HALYARD (844-425-9273)** or email MyHalyardOrders@owens-minor.com