

Support Bandages

Tubular Elastic Support Bandages

- Tubular elastic, cotton/spandex bandage for support and compression
- Easy to apply and stays in place without pins, clips or tape
- Delivers consistent support with a single application
- Convenient alternative to flat wrap bandages

Product Code	Product Description	Quantity
TSB5013	Bandage Support Natural 13in x 11yd Large Trunk Size L	1 Each/Case
TSB5150	Bandage Support Natural 1-1/2in x 11yd Infant Feet And Arms Size A	1 Each/Case
TSB5250	Bandage Support Natural 2-1/2in x 11yd Small Hands And Arms Size B	1 Each/Case
TSB5275	Bandage Support Natural 2-3/4in x 11yd Med Arm Small Ankle Size C	1 Each/Case
TSB5300	Bandage Support Natural 3in x 11yd Large Arms Med Ankles Size D	1 Each/Case
TSB5350	Bandage Support Natural 3-1/2in x 11yd Lg Ankle Med Knee Size E	1 Each/Case
TSB5400	Bandage Support Natural 4in x 11yd Large Knee Med Thigh Size F	1 Each/Case
TSB5450	Bandage Support Natural 4-1/2in x 11yd Large Thigh Size G	1 Each/Case
TSB5550	Bandage Support Natural 5-1/2in x 11yd Large Knee XL Size I	1 Each/Case
TSB5675	Bandage Support Natural 6-3/4in x 11yd Small Trunk Size J	1 Each/Case
TSB5825	Bandage Support Natural 8-1/4in x 11yd Medium Trunk Size K	1 Each/Case
TSB7250	Bandage Support Beige 2-1/2in x 11yd Small Hands And Arms Size B	1 Each/Case
TSB7275	Bandage Support Beige 2-3/4in x 11yd Med Arm Small Ankle Size C	1 Each/Case
TSB7300	Bandage Support Beige 3in x 11yd Large Arms Med Ankles Size D	1 Each/Case
TSB7350	Bandage Support Beige 3-1/2in x 11yd Lg Ankle Med Knee Size E	1 Each/Case
TSB7400	Bandage Support Beige 4in x 11yd Large Knee Med Thigh Size F	1 Each/Case
TSB7450	Bandage Support Beige 4-1/2in x 11yd Large Thigh Size G	1 Each/Case
TSB9250	Bandage Support Natural 2-1/2in x 36in Small Hands And Arms Size B	12 Each/Case
TSB9275	Bandage Support Natural 2-3/4in x 36in Med Arm Small Ankle Size C	12 Each/Case
TSB9300	Bandage Support Natural 3in x 36in Large Arms Med Ankles Size D	12 Each/Case
TSB9350	Bandage Support Natural 3-1/2in x 36in Leg Ankle Med Knee Size E	12 Each/Case
TSB9400	Bandage Support Natural 4in x 36in Large Knee Medium Thigh Size F	12 Each/Case



MediChoice® Tubular Elastic Support Bandage Compression Sizing Chart

Tubular Elastic Support Bandage Compression Sizing Chart				Low Pressure 5–10 mmHG		Medium Pressure 10–20 mmHG		High Pressure 20–30 mmHG	
Catalog#	Color	Size Code	Size	Centimeters	Inches	Centimeters	Inches	Centimeters	Inches
TSB5150	Natural	Size A	1.50"	19-22 cm	7.5-8.7 in	22-24 cm	8.7-9.4 in	24-27 cm	9.4-10.6 in
TSB5250	Natural	Size B	2.50"	22-24 cm	8.7-9.4 in	24-27 cm	9.4-10.6 in	27-30 cm	10.6-11.8 in
TSB5275	Natural	Size C	2.75"	24-27 cm	9.4-10.6 in	27-30 cm	10.6-11.8 in	30-35 cm	11.8-13.8 in
TSB5300	Natural	Size D	3.00"	27-30 cm	10.6-11.8 in	30-35 cm	11.8-13.8 in	35-40 cm	13.8-17.3 in
TSB5350	Natural	Size E	3.50"	30-35 cm	11.8-13.8 in	35-44 cm	13.8-17.3 in	40-48 cm	17.3-18.9 in
TSB5400	Natural	Size F	4.00"	35-44 cm	13.8-17.3 in	44-58 cm	17.3-22.8 in	58-64 cm	22.8-25.2 in
TSB5450	Natural	Size G	4.50"	44-58 cm	17.3-22.8 in	58-64 cm	22.8-25.2 in	64-72 cm	25.2-28.3 in
TSB5550	Natural	Size I	5.50"	44-58 cm	17.3-22.8 in	58-64 cm	22.8-25.2 in	64-72 cm	25.2-28.3 in
TSB5675	Natural	Size J	6.75"	64-72 cm	25.2-28.3 in	72-86 cm	28.3-33.8 in	86-96 cm	33.8-37.8 in
TSB5825	Natural	Size K	8.25"	72-86 cm	28.3-33.8 in	86-96 cm	33.8-37.8 in	96-108 cm	37.8-42.5 in
TSB5013	Natural	Size L	8.25"	72-86 cm	28.3-33.8 in	86-96 cm	33.8-37.8 in	96-108 cm	37.8-42.5 in
TSB7250	Beige	Size B	2.50"	22-24 cm	8.7-9.4 in	24-27 cm	9.4-10.6 in	27-30 cm	10.6-11.8 in
TSB9250	Beige	Size B	2.50"	22-24 cm	8.7-9.4 in	24-27 cm	9.4-10.6 in	27-30 cm	10.6-11.8 in
TSB7275	Beige	Size C	2.75"	24-27 cm	9.4-10.6 in	27-30 cm	10.6-11.8 in	30-35 cm	11.8-13.8 in
TSB9275	Beige	Size C	2.75"	24-27 cm	9.4-10.6 in	27-30 cm	10.6-11.8 in	30-35 cm	11.8-13.8 in
TSB7300	Beige	Size D	3.00"	27-30 cm	10.6-11.8 in	30-35 cm	11.8-13.8 in	35-40 cm	13.8-17.3 in
TSB9300	Beige	Size D	3.00"	27-30 cm	10.6-11.8 in	30-35 cm	11.8-13.8 in	35-40 cm	13.8-17.3 in
TSB7350	Beige	Size E	3.50"	30-35 cm	11.8-13.8 in	35-44 cm	13.8-17.3 in	40-48 cm	17.3-18.9 in
TSB9350	Beige	Size E	3.50"	30-35 cm	11.8-13.8 in	35-44 cm	13.8-17.3 in	40-48 cm	17.3-18.9 in
TSB7400	Beige	Size F	4.00"	35-44 cm	13.8-17.3 in	44-58 cm	17.3-22.8 in	58-64 cm	22.8-25.2 in
TSB9400	Beige	Size F	4.00"	35-44 cm	13.8-17.3 in	44-58 cm	17.3-22.8 in	58-64 cm	22.8-25.2 in
TSB7450	Beige	Size G	4.50"	44-58 cm	17.3-22.8 in	58-64 cm	22.8-25.2 in	64-72 cm	25.2-28.3 in

