


3-Layer Compression Bandage System with Trak™ Bandage

- 3-Layer system for optimum compression and absorbency in the treatment of venous leg ulcers of varying stages and other related conditions
 - Layer 1: Non-woven absorbent padding for comfort and protection of bony joints
 - Layer 2: Patented Trak™ compression bandage provides a calculated guide for clinicians to deliver consistent compression with every application
 - Layer 3: Self-adhering cohesive bandage secures system together
- Pressure: 35–40 mm Hg
- Three individually packaged layers in a resealable pouch
- Adhesive tape strips and detailed instruction guide included

NOT MADE WITH
natural rubber latex

diverse
supplier

Sterile Premium Dual Self-Closure		
Catalog #	Description	Quantity
BND3040	3-Layer System With Trak Guide™	8 EA/CA



3-Layer Reduced Compression Bandage System

- 3-Layer system for optimum compression and absorbency in the treatment of venous leg ulcers of varying stages
 - Layer 1: Non-woven absorbent padding for comfort and protection of bony joints
 - Layer 2: Breathable cotton/poly crepe conforming bandage to anchor and smooth padded layer while absorbing exudate
 - Layer 3: Self-adhering cohesive bandage secures system together
- Pressure: 17-25 mm/Hg
- Three individually packaged layers in a resealable pouch
- Adhesive tape strips and detailed instruction guide included

NOT MADE WITH
natural rubber latex

diverse
supplier

Ordering Information		
Catalog #	Description	Quantity
BND3000	3-Layer System	8 EA/CA



4-Layer Compression Bandage System

- Four Layer Bandage Kit is designed for optimum compression and absorbency in the treatment of venous leg ulcers and related conditions
 - Layer 1: Non-woven absorbent padding for comfort and protection of bony joints
 - Layer 2: Breathable cotton/poly crepe conforming bandage to anchor and smooth padded layer while absorbing exudate
 - Layer 3: Light compression bandage conforms to any contour while providing half of the compression of the system
 - Layer 4: Self-adhering cohesive bandage secures system together
- Pressure: 35-40 mm/Hg
- Four individually packaged layers in a resealable pouch
- Adhesive tape strips and detailed instruction guide included

NOT MADE WITH
natural rubber latex

diverse
supplier

Ordering Information		
Catalog #	Description	Quantity
BND9001	4-Layer System	8 EA/CA

