




PRIMO Glide & Guide Turn Sheet

PRIMO Glide & Guide Turn Sheet allows for easy gliding across the bed while reducing staff injuries.

PRIMO'S Glide & Guide Turn Sheet **reduces friction by up to 67%** significantly reducing the force needed to position patients and offload the sacrum, reposition patients with ease, and simplify lateral transfers.

Minimize risk of friction and skin shear with high quality patient side fabric, micro-suede



Allows staff members to move the patient in bed with top and bottom handles




Slides with ease, reducing effort needed to reposition patient and decreasing back and shoulder injuries



Help reduce the sliding/bunching of the sheet beneath the patient with top strap and side strap options



Ergonomically turn patient with ease with roll strap options



Allows for convenient gripping with x-box stitching, and eight 1½" handles on each side of the sheet



Available in a variety of sizes offering various widths and lengths for bariatric and standard beds

Turn Sheet Weight Limits
Standard: < 565 lbs
Bariatric: < 1,000 lbs

PRIMO Glide & Guide Body Wedge

PRIMO Glide & Guide Body Wedges reliably offload the sacrum while maintaining a consistent 30° patient position.

Non-slip bottom keeps wedge in place

Fluid resistant cover allows for a completely wipeable surface

Handles are convenient for gripping to remove and reposition patient



Offloads sacrum with a 30° side-lying angle position

PRIMO Glide & Guide Body Wedges

Wedges are available with a classic angled edge or rounded edge. Rounded body wedge decreases pressure on the back for **increased patient comfort**.

PRIMO Underpad



Patent Pending

PRIMO Glide & Guide Turn System



One Sheet + Two Wedges + Underpad = One Bundle

Bundles (Kits) All bundles include one sheet and two body wedges							
Item #	Product Description	Size	Case QTY	Sheet Included	Underpad Included	Wedge Set Included	Bundle Description
GGTS13	Glide & Guide Turn System Angled	Standard	10 Kits	GGGS3660TBH	UPA3657	BPW11F	Turn Sys (1) 36x60" Sheet, (2) Wdg, (5) Upads
GGTS14-5	Glide & Guide Turn System Angled	Bariatric	5 Kits	GGGS5670TBH	UPA4057	BPW11F	Turn Sys (1) 56x70" Sheet, (2) Wdg, (5) Upads
GGTS09-5	Glide & Guide Turn System Rounded	Standard	5 Kits	GGGS3660	UPA3657	BPW11R	Turn Sys (1) 36x60" Sheet, (2) Wdg, (5) Upads
GGTS09B-5	Glide & Guide Turn System Rounded	Bariatric	5 Kits	GGGS5660	UPA4057	BPW11R	Turn Sys (1) 56x60" Sheet, (2) Wdg, (5) Upads

Turn Sheets

Item #	Description	Color	W x L	Size	Case QTY
GGS3660TBH	Turn Sheet - Top/Bottom Handles	Purple/Gray	36 x 60	Standard	24 EA
GGS5670TBH	Turn Sheet - Top/Bottom Handles	Navy/Gray	56 x 70	Bariatric	24 EA
GGS3660	Standard Turn Sheet	Purple/Gray	36 X 60	Standard	24 EA
GGS5660	Bariatric Turn Sheet	Green/Gray	56 X 60	Bariatric	24 EA

Eaches



Angled



Rounded

Eaches - Wipe Down Body Wedges

Item #	Edge	Length	Case QTY
BPW04FL	Angled	19" Long	10 each
BPW04LR	Rounded	19" Long	10 each

Sets



Sets - Wipe Down Body Wedges

Item #	Edge	Length	Case QTY
BPW11F	Angled	2 each - 15" Long	12 sets
BPW11R	Rounded	2 each - 15" Long	10 sets

Underpads

Underpads

Item #	Description	Length	Case QTY
UPA3657	Super Absorbency Underpad, Standard	36 x 57	25EA - 5EA/BG 5BG/CS
UPA4057	Super Absorbency Underpad, Bariatric	40 x 57	25EA - 5EA/BG 5BG/CS