

FEMORAL ANGIOGRAPHY DRAPES


DRAPES DESIGNED TO PROTECT YOU AND YOUR PATIENTS IN THE CATH LAB

In high-fluid diagnostic and interventional procedures performed in the Catheterization and Interventional Radiology labs, effective patient access, fluid control and full coverage are critical for the protection of both patients and staff. Now you have more choice than ever with five HALYARD* Drapes that deliver performance, protection and convenience to meet the special needs of these procedures.

THE RIGHT DRAPE FOR EVERY PROCEDURE

As a recognized leader in the development of high-performance medical fabrics, Halyard has rights to over 600 patents for advanced nonwovens materials. Like all of our drapes, our newest additions provide the critical balance of properties your demanding work requires:¹

- CONTROL* Plus reinforcement fabric provides the highest rating for flame resistance¹ and resistance to laser ignition².
- Low-lint³ and abrasion-resistant fabric⁴.
- Fabric resistant to tearing⁵ and strikethrough⁶.
- Surgical Drapes with CONTROL* Plus reinforcement fabric meets AAMI Level 4 liquid barrier performance test requirements⁶.



INTRODUCING FIVE FEMORAL ANGIOGRAPHY DRAPES


FEMORAL ANGIOGRAPHY/IJ ACCESS COMBINATION DRAPE

The newest addition to our line of femoral angiography drapes gives you a choice of four fenestrations, allowing for femoral and/or IJ access on either side of the patient's body, with clear film windows.

- Two 3.5" adhesive circular fenestrations
- Two 7" x 7" trapezoidal fenestrations
- Absorbent CONTROL* Plus reinforcement fabric
- Clear side windows to allow visibility of table controls

Code	Description
29225	Femoral Angiography/IJ Access Combination Drape
29226NS	Femoral Angiography/IJ Access Combination Drape

All codes in bold indicate a sterile product
Non-sterile product available for custom trays.




FEMORAL ANGIOGRAPHY DRAPE WITH SIDE & BOTTOM WINDOWS

This new femoral angiography drape features the extra-long coverage preferred for cath lab procedures, plus clear windows along the sides and the end of the drape to allow maximum visibility for controls.

- Two 3.5" adhesive circular fenestrations
- Clear film side and bottom windows
- Absorbent CONTROL* PLUS reinforcement fabric

Code	Description
29227	Femoral Angiography Drape with Side & Bottom Windows
29228NS	Femoral Angiography Drape with Side & Bottom Windows

All codes in bold indicate a sterile product
Non-sterile product available for custom trays.




PEDIATRIC FEMORAL ANGIOGRAPHY DRAPE

Sized especially for smaller patients, our Pediatric Femoral Angiography drape includes the same features as our adult/ full-size drapes, including a choice of four fenestrations for optimal access.

- Two 1.75" adhesive circular fenestrations
- Two 3.5" adhesive circular fenestrations
- Dual clear film poly windows
- Absorbent CONTROL* PLUS reinforcement fabric

Code	Description
29232	Pediatric Femoral Angiography Drape
29233NS	Pediatric Femoral Angiography Drape

All codes in bold indicate a sterile product
 Non-sterile product available for custom trays.




PACEMAKER DRAPE WITH FENESTRATED INCISE

Our newest drape for pacemaker procedures features two clear incise windows with adhesive fenestrations to aid correct placement.

- Two 3" x 5" fenestrated incise areas
- Dual clear film poly windows
- 10" Fabric trough (made with absorbent CONTROL* Plus reinforcement fabric)
- Fabric tube holders
- Absorbent CONTROL* PLUS reinforcement fabric

Code	Description
29229	Pacemaker Drape with Fenestrated Incise
29231NS	Pacemaker Drape with Fenestrated Incise

All codes in bold indicate a sterile product
 Non-sterile product available for custom trays.



FEMORAL ANGIOGRAPHY DRAPES

5 Digit Code	Type	Description	Case Count
29225	Sterile	Femoral Angiography/IJ Access Combination Drape	15
29226NS	Non-sterile	Femoral Angiography/IJ Access Combination Drape	15
29227	Sterile	Femoral Angiography Drape with Side and Bottom Windows	13
29228NS	Non-sterile	Femoral Angiography Drape with Side and Bottom Windows	14
29232	Sterile	Pediatric Femoral Angiography Drape	12
29233NS	Non-sterile	Pediatric Femoral Angiography Drape	15
29229	Sterile	Pacemaker Drape with Fenestrated Incise	12
29231NS	Non-sterile	Pacemaker Drape with Fenestrated Incise	12
89583	Sterile	HALYARD BASICS* Radial Femoral Angiography Drape with Dual Windows	12
79583	Non-sterile	HALYARD BASICS* Radial Femoral Angiography Drape with Dual Windows	14

All codes in **bold** indicate a sterile product. Non-sterile product available for custom trays.




KNOWLEDGE NETWORK* Accredited Education
In-Service Training and Technical Support
Credentialed Sales representatives
Tools & Best Practices
Clinical Research
Commitment to Excellence

For more information, please visit:
www.halyardhealth.com

Call 1-844-HALYARD (1-844-425-9273)
in the United States and Canada.



*Registered Trademark or Trademark of Halyard Health, Inc.
© 2015 HYH. All rights reserved. C15553



HALYARD BASICS* RADIAL FEMORAL ANGIOGRAPHY DRAPE

Our drape for radial access gives you a choice of four fenestrations, allowing for radial or femoral access on either side of the patient's body, with adhesive fenestrations

- Two 4.72" adhesive circular fenestrations
- Two oval 4.72" x 2.7" radial fenestrations
- Dual clear film

Code	Description
89583	HALYARD BASICS* Radial Femoral Angiography Drape with Dual Windows
79583	HALYARD BASICS* Radial Femoral Angiography Drape with Dual Windows

All codes in **bold** indicate a sterile product
Non-sterile product available for custom trays.

Ask about HALYARD* C-Arm and Banded Bag Equipment Covers!

- 1 Flammability Standard: 16FR 1610: 2003-Class 1
 - 2 ISO 11810: 2005 Part One
 - 3 Test method ISO 9073-10; Criteria <4.0 Log10 (lint count).
 - 4 Test method ASTM D4970:2007 Criteria ≥ 3
 - 5 Test method DSS87-08:2008
 - 6 Test method ASTM F1670: 2008
- + Claims do not apply to HALYARD BASICS* Radial Femoral Angiography Drape