



MULTI-ENZYME DETERGENT

THE INSTRUMENT CLEANER YOU NEED FOR BETTER REPROCESSING OUTCOMES

Proper cleaning of medical devices is the critical step in decontamination, prior to disinfection or sterilization. HALYARD* Multi-Enzyme Detergent rapidly removes contaminants from reusable medical devices used throughout your facility.

The synergistic formulation of enzymes plus detergent works quickly, over a range of water temperatures. It penetrates and breaks down the most common contaminants found on your devices, including protein, fat and carbohydrate-based soils. HALYARD* Multi-Enzyme Detergent helps assure your devices are clean before being disinfected or sterilized, which can result in better reprocessing outcomes and reduced risk of cross-contamination from instruments.

HALYARD* Multi-Enzyme Detergent provides benefits that can help make your facility's cleaning process easier and more complete:

- Effective: rapidly removes all types of soils
- Compatible: can be used with a wide range of devices
- Low-Foaming: appropriate for manual and small washer use
- Easy-To-Use: pump dispenses consistent 1 ounce dose
- User-Friendly: pleasant fragrance in an ergonomically designed package
- Process-Friendly: easily integrates into your current practices
- Environment-Friendly: biodegradable and phosphate-free

HALYARD* Multi-Enzyme Detergent is just one of the clinical solutions you can depend on to meet the demands of your fast-paced world.

For more information about HALYARD* Multi-Enzyme Detergent, ask your Halyard Sales Representative, authorized Halyard Distributor Representative or visit our website at www.halyardhealth.com



HALYARD* MULTI-ENZYME DETERGENT



HALYARD* MULTI-ENZYME DETERGENT

Stock#	Description	Case Information
65201 1	Gallon Bottle	4 (1 gallon) bottles w/1 dispensing pump

Also available: 24 oz. HALYARD* Pre-Cleaning Detergent Spray

PPE PRODUCTS

- Fluid-Resistant Facial and Eye Protection
- Impervious Cover Gowns
- Nitrile Exam Gloves
- Bouffant Caps
- Shoe Covers

ENDOSCOPY PRODUCTS

- Endoscope Cleaning Brushes
- Endoscope Cleaning Sponges

EDUCATIONAL PROGRAMS

"Coming Clean: An Essential for Reusable Medical Devices" An educational video presentation, CE-accredited for both Nurses and CS/SPD Technicians. Ask your Halyard Sales Representative to schedule a presentation at your facility.



KNOWLEDGE NETWORK* Accredited Education
In-Service Training and Technical Support
Credentialed Sales representatives
Tools & Best Practices
Clinical Research
Commitment to Excellence

For more information, please visit:
www.halyardhealth.com

Call 1-844-HALYARD (1-844-425-9273)
in the United States and Canada.



*Registered Trademark or Trademark of Halyard Health, Inc.
© 2015 HYH. All rights reserved. C15562