

**HALYARD\***

Two-Tier Back Table Covers



# BETTER COVERAGE. BETTER VALUE.

**TWO LEVELS OF PROTECTION:** A tough, 5.5 mil poly film that resists holes and tears and a CONTROL\* absorbent reinforcement layer.

- The heavy-duty film layer is designed to stand up to your heaviest trays and instruments.
- The CONTROL\* layer gives you a level of absorbency that may reduce your need for added OR towels on the back table. Along with saving time and money, using fewer towels helps keep lint out of the OR, which decreases the risk of wound contamination.<sup>1</sup>



## BETTER BACK TABLE COVERAGE. BETTER UP-FRONT VALUE.

### A NEW FOLD: Supports aseptic placement.


Completely film-reinforced and fluid-resistant throughout; our strong fabric is flexible, drapable, and easy to aseptically position on your two-tier back table.

### THREE SIZES: Fits two-tier tables 48", 60" and 72".

Our one-piece construction features a proprietary<sup>2</sup> design that makes it simple to fit two-tier back tables currently on the market. Extra reinforced material at the table edges protects against tearing. The proprietary non-slip feature helps keep the cover in place while in use.

### HALYARD\* TWO TIER PADDED HEAVY DUTY BACK TABLE COVERS

| Product Codes |             |                                   | Case Count |             |
|---------------|-------------|-----------------------------------|------------|-------------|
| Sterile       | Non-Sterile | Cover Fits:                       | Sterile    | Non-Sterile |
| 37496         | 37495       | 48" / 121 cm Two-Tier Back Tables | 10         | 12          |
| 37498         | 37497       | 60" / 152 cm Two-Tier Back Tables | 10         | 12          |
| 37501         | 37499       | 72" / 182 cm Two-Tier Back Tables | 10         | 12          |



### CUSTOM TRAYS?

Ask for HALYARD\* 2-TIER BACK TABLE COVERS in your next CPT!



KNOWLEDGE NETWORK\* Accredited Education  
In-Service Training and Technical Support  
Credentialed Sales Representatives  
Tools & Best Practices  
Clinical Research  
Commitment to Excellence



For more information, please visit:

[www.halyardhealth.com](http://www.halyardhealth.com)

Call 1-844-HALYARD (1-844-425-9273)  
in the United States and Canada.

[www.halyardhealth.com](http://www.halyardhealth.com)



1. Janoff K. Foreign Body Reactions Secondary to cellulose lint fibers. American Journal of Surgery, vol: 147, NOS, pp 598-600, 1984.  
2. Patent pending.

\*Registered Trademark or Trademark of Halyard Health, Inc. or its affiliates. © 2014 HYH. All rights reserved.  
SS63 C14181