



# TECHNICAL BULLETIN

## Clarification of Product Storage Requirements ON-Q\* SilverSoaker\* Catheters

The ON-Q\* SilverSoaker\* Catheters shipping box label and sterile package label contain the following symbol regarding product storage requirements:



This symbol is defined in the International Guidance Document ISO15223-1:2012 (Symbols to be Used with Medical Device Labels, Labeling and Information to be Supplied) and has the following meaning:

SYMBOL	SYMBOL TITLE	MEANING OF SYMBOL
	Keep away from sunlight	Indicates a medical device that needs protection from light sources.

Although the title of the symbol is "Keep away from sunlight," the description of the symbol "Indicates a medical device that needs protection from light sources."

Halyard recommends the following practices for protecting the product from light sources:

- Store ON-Q\* SilverSoaker\* catheters in an area away from all direct light sources. Storage within the shipping box (i.e., corrugate box) is best.

In addition:

- Ensure that products are used prior to the expiration date shown on the package label.
- Refer to the ON-Q\* SilverSoaker\* Instructions for Use (IFU) for complete instructions, cautions and warnings.

Materials Managers are encouraged to inform all users of ON-Q\* SilverSoaker\* Catheters within your facility about the appropriate storage requirements.

There are inherent risks in all medical devices. Please refer to the product labeling for **Indications, Cautions, Warnings and Contraindications**. Failure to follow the product labeling could directly impact patient safety. Physician is responsible for prescribing and administering medications per instructions provided by the drug manufacturer. Refer to [www.halyardhealth.com](http://www.halyardhealth.com) for additional product safety **Technical Bulletins**.

Please contact the Clinical Services Department at 800-444-2728 or 949-923-2400 if you have any questions regarding this information.