

# YOU'RE FIGHTING CANCER.

## Wear the Battle Gear That Protects Against 52 Chemotherapy Drugs!



### It's 52 vs. You.

PURPLE NITRILE-XTRA\* Gloves and HALYARD\* Procedure Gowns for Chemotherapy Use have been tested for **resistance to 52 chemotherapy drugs.**<sup>†</sup>

Arsenic Trioxide (1 mg/ml)	Cyclophosphamide (20 mg/ml)	Fludarabine (25 mg/ml)	Mitomycin (0.5 mg/ml)	Temsirolimus (25 mg/ml)
Azacitidine (Vidaza) (25 mg/ml)	Cytarabine HCl (100 mg/ml)	Fluorouracil (50 mg/ml)	Mitoxantrone (2 mg/ml)	Trastuzumab (21 mg/ml)
Bendamustine (5 mg/ml)	Cytovene (10 mg/ml)	Fulvestrant (50 mg/ml)	Oxaliplatin (2 mg/ml)	ThioTEPA (10 mg/ml)
Bleomycin Sulfate (15 mg/ml)	Dacarbazine (10 mg/ml)	Gemcitabine (38 mg/ml)	Paclitaxel (6 mg/ml)	Topotecan HCl (1 mg/ml)
Bortezomib (Velcade) (1 mg/ml)	Daunorubicin HCl (5 mg/ml)	Idarubicin (1 mg/ml)	Paraplatin (10 mg/ml)	Triclosan (1 mg/ml)
Busulfan (6 mg/ml)	Decitabine (5 mg/ml)	Ifosfamide (50 mg/ml)	Pemetrexed (25 mg/ml)	Trisenox (0.1 mg/ml)
Carboplatin (10 mg/ml)	Docetaxel (10 mg/ml)	Irinotecan (20 mg/ml)	Pertuzumab (30 mg/ml)	Vincristine Sulfate (1 mg/ml)
Carfilzomib (2 mg/ml)	Doxorubicin HCl (2 mg/ml)	Mechlorethamine HCl (1 mg/ml)	Raltitrexed (0.5 mg/ml)	Vinblastine (1 mg/ml)
Carmustine (3.3 mg/ml) <sup>†</sup>	Ellence (2 mg/ml)	Melphalan (5 mg/ml)	Retrovir (10 mg/ml)	Vinorelbine (10 mg/ml)
Cetuximab (Erbix) (2 mg/ml)	Eribulin Mesylate (0.5 mg/ml)	Methotrexate (25 mg/ml)	Rituximab (10 mg/ml)	Zoledronic Acid (0.8 mg/ml)
Cisplatin (1 mg/ml)	Etoposide (20 mg/ml)			

<sup>†</sup> Testing measured no breakthrough at the Standard Breakthrough Rate of 0.1ug/cm<sup>2</sup>/minute, up to 240 minutes for gloves and 480 minutes for gowns, except for carmustine. Carmustine breakthrough occurred at 80.4 minutes on gloves and 33.3 minutes on gowns.

# There is no safe level of exposure<sup>1</sup> to chemotherapy drugs for healthcare workers.

Is your team at risk? New, stricter guidelines — **USP <800>** — apply to *all healthcare personnel* who handle hazardous drugs.<sup>2</sup>

Don't wait: get full protection now with **HALYARD\* Chemo Gear!**

## **PURPLE NITRILE-XTRA\*** Exam Gloves

Code	Description	Case Count
55090	X-Small, non-sterile, ambidextrous	500
50601	Small, non-sterile, ambidextrous	500
50602	Medium, non-sterile, ambidextrous	500
50603	Large, non-sterile, ambidextrous	500
50604	X-Large, non-sterile, ambidextrous	500

More available at [halyardhealth.com](http://halyardhealth.com).

## **NEW PINK UNDERGUARD\*** Nitrile Exam Gloves

Code	Description	Case Count
47452	X-Small, non-sterile, ambidextrous	1000
47453	Small, non-sterile, ambidextrous	1000
47454	Medium, non-sterile, ambidextrous	1000
47455	Large, non-sterile, ambidextrous	1000
47456	X-Large, non-sterile, ambidextrous	1000

More available at [halyardhealth.com](http://halyardhealth.com).

## **HALYARD\* FACE SHIELDS**

Code	Description	Case Count
41204	Full length clear visor, 9.5"	40
41205	Three-quarter length clear visor, 7.5"	40

## **NEW CHEMO360\***

Procedure Gowns Tested for Chemotherapy Use

Code	Description	Case Count
47346	Universal, Closed back, Tape closure, Poly-coated	100
47347	XXL, Closed back, Tape closure, Poly-coated	100

## **HALYARD\* PROCEDURE GOWNS** Tested for Chemotherapy Use

Code	Description	Case Count
69606	Universal, Open back, Over-the-head, Poly-coated	100
37284	XX Large, Over-the-head, Poly-coated	100

## **FLUIDSHIELD\*** N95 RESPIRATOR MASKS

Code	Description	Case Count
46828	Small, individually bagged	210
46728	Regular, individually bagged	210
46827	Small	210
46727	Regular	210
46867	Small, with SAFETY SEAL*	210
46767	Regular, with SAFETY SEAL*	210

More available at [halyardhealth.com](http://halyardhealth.com).



Get the facts. Visit  
[www.GearUpForTheFight.com](http://www.GearUpForTheFight.com)

Call 1-844-HALYARD (1-844-425-9273)  
in the United States and Canada.

[www.halyardhealth.com](http://www.halyardhealth.com)



1. Centers for Disease Control and Prevention, Hazardous Drug Exposures in Health Care, [www.cdc.gov/-niosh/topics/hazdrug/default.html](http://www.cdc.gov/-niosh/topics/hazdrug/default.html)  
2. USP General Chapter <800> Hazardous Drugs – Handling in Healthcare Settings, United States Pharmacopoeial Convention, USP 40–NF 35, Second Supplement (2017), 1. [www.usp.org](http://www.usp.org)

\* Registered Trademark or Trademark of O&M Halyard or its affiliates. © 2018. All rights reserved.  
\* PURPLE, PURPLE NITRILE and PURPLE NITRILE X-TRA are Registered Trademarks or Trademarks of O&M Halyard or its affiliates.  
The COLOR PURPLE is a Registered Trademark or Trademark O&M Halyard or its affiliates.  
C182481-E