

# Bandages and Rolls

## Conforming Stretch Gauze Bandages

- Soft, one-ply polyester knit weave
- Breathable, absorbent gauze stretches and clings to itself while providing soft cushioning
- Finished edges minimize linting and unraveling
- Ideal secondary dressing gently conforms to body contours without restricting movement
- Can be applied stretched or left unstretched; length is measured unstretched
- Available in a range of widths and lengths as well as soft packs and bulk bags



Product Code	Product Description	Quantity
CFBN1001	Conforming Bandage 1in x 75in Stretch Gauze Sterile	1 Each, 12 Each/Box, 8 Boxes/Case, 96 Each/Case
CFBN1002	Conforming Bandage 2in x 75in Stretch Gauze Sterile	1 Each, 12 Each/Box, 8 Boxes/Case, 96 Each/Case
CFBN1003	Conforming Bandage 3in x 75in Stretch Gauze Sterile	1 Each, 12 Each/Box, 8 Boxes/Case, 96 Each/Case
CFBN1004	Conforming Bandage 4in x 75in Stretch Gauze Sterile	1 Each, 12 Each/Box, 8 Boxes/Case, 96 Each/Case
CFBN1006	Conforming Bandage 6in x 82in Stretch Gauze Sterile	1 Each, 12 Each/Box, 4 Boxes/Case, 48 Each/Case
CFBN2001	Conforming Bandage 1in x 75in Stretch Gauze Non Sterile	24 Each/Bag, 4 Bags/Case
CFBN2002	Conforming Bandage 2in x 75in Stretch Gauze Non Sterile	12 Each/Bag, 8 Bags/Case
CFBN2003	Conforming Bandage 3in x 75in Stretch Gauze Non Sterile	12 Each/Bag, 8 Bags/Case
CFBN2004	Conforming Bandage 4in x 75in Stretch Gauze Non Sterile	12 Each/Bag, 8 Bags/Case
CFBN2006	Conforming Bandage 6in x 82in Stretch Gauze Non Sterile	6 Each/Bag, 8 Bags/Case
CFBN2007	Conforming Bandage 2in x 36in Stretch Gauze Non Sterile	100 Each/Bag, 8 Bags/Case

## Fluff Gauze Bandage Rolls

- 100% cotton gauze, washed and fluff-dried
- Open, crimped weave for fluid absorption, aeration and comfort
- Bulk, loft and light stretch to conform to body contours and provide cushioning
- Finished edges minimize linting and unraveling
- Used as a primary absorbent dressing layer or as a secondary dressing to secure and protect primary dressings
- Ideal for bandaging head, limbs and difficult-to-dress sites
- Available in 3-ply, 6-ply, and 8-ply
- Available in a range of sizes



Product Code	Product Description	Quantity
GZBN1001	Gauze Bandage Roll 2in x 4yd Small 3-Ply Sterile	1 Each, 96 Each/Case
GZBN1002	Gauze Bandage Roll 3in x 4yd Medium 3-Ply Sterile	1 Each, 96 Each/Case
GZBN1003	Gauze Bandage Roll 4in x 4yd Large 3-Ply Sterile	1 Each, 96 Each/Case
GZBN1004	Gauze Bandage Roll 6in x 4yd X-Large 3-Ply Sterile	1 Each, 48 Each/Case
GZBN2001	Gauze Bandage Roll 2in x 4yd Small 3-Ply Non Sterile	12 Each/Bag, 8 Bags/Case
GZBN2002	Gauze Bandage Roll 3in x 4yd Medium 3-Ply Non Sterile	12 Each/Bag, 8 Bags/Case
GZBN2003	Gauze Bandage Roll 4in x 4yd Large 3-Ply Non Sterile	12 Each/Bag, 8 Bags/Case
GZBN2004	Gauze Bandage Roll 6in x 4yd X-Large 3-Ply Non Sterile	6 Each/Bag, 8 Bags/Case
GZBN3001	Gauze Bandage Roll 2-1/4in x 3yd Small 6-Ply Sterile	1 Each, 96 Each/Case
GZBN3002	Gauze Bandage Roll 3-13/32in x 3-39/64yd Medium 6-Ply Sterile	1 Each, 96 Each/Case
GZBN3003	Gauze Bandage Roll 4-1/2in x 4-7/64yd Large 6-Ply Sterile	1 Each, 100 Each/Case

Product Code	Product Description	Quantity
GZBN3004	Gauze Bandage Roll 4-1/2in x 4-4/64yd Large 6-Ply Sterile	1 Each/Tray, 60 Trays/Case
GZBN4001	Gauze Bandage Roll 2-1/4in x 3yd Small 6-Ply Non Sterile	12 Each/Bag, 8 Bags/Case
GZBN4002	Gauze Bandage Roll 3-13/32in x 3-39/64yd Medium 6-Ply Non Sterile	96 Each/Case
GZBN4003	Gauze Bandage Roll 4-1/2in x 4-7/64yd Large 6-Ply Non Sterile	100 Each/Case
GZBN4004	Gauze Bandage Roll 4-1/2in x 4-7/64yd Large 6-Ply Bag Non Sterile	12 Each/Bag, 4 Bags/Case
GZBN5001	Gauze Bandage Roll 4-1/2in x 3-7/64yd Large 8-Ply Sterile	1 Each, 100 Each/Case