

Non-Adherent Dry Dressing

Adhesive Island Dressings

- General purpose bordered gauze dressing
- Easy to use all-in-one dressing is conformable and a gentle alternative to gauze and tape
- Non-adherent contact layer with absorptive pad
- Soft, absorbent pad allows light to moderate fluid absorption without adhering to the wound or disrupting healing tissue
- Non-woven, breathable and non-irritating adhesive border with rounded edges secures the dressing in place and protects the wound
- Split release liner for easy, one-step application
- Ideal for incision sites and other wounds as a primary or secondary dressing
- Available in a range of sizes
- Available in standard as well as plus option for increased fluid absorption



Product Code	Product Description	Quantity
IDSG1001	Island Dressing 2in x 3in Standard Adhesive Sterile	1 Each, 50 Each/Box, 8 Boxes/Case, 400 Each/Case
IDSG1002	Island Dressing 4in x 4in Standard Adhesive Sterile	1 Each, 25 Each/Box, 8 Boxes/Case, 200 Each/Case
IDSG1003	Island Dressing 4in x 5in Standard Adhesive Sterile	1 Each, 25 Each/Box, 8 Boxes/Case, 200 Each/Case
IDSG1004	Island Dressing 4in x 8in Standard Adhesive Sterile	1 Each, 25 Each/Box, 4 Boxes/Case, 100 Each/Case
IDSG1005	Island Dressing 4in x 10in Standard Adhesive Sterile	1 Each, 25 Each/Box, 4 Boxes/Case, 100 Each/Case
IDSG1006	Island Dressing 4in x 14in Standard Adhesive Sterile	1 Each, 25 Each/Box, 2 Boxes/Case, 50 Each/Case
IDSG1007	Island Dressing 6in x 6in Standard Adhesive Sterile	1 Each, 25 Each/Box, 4 Boxes/Case; 100 Each/Case
IDSG2001	Island Dressing 4in x 6in Plus Adhesive Sterile	1 Each, 25 Each/Box, 4 Boxes/Case, 100 Each/Case
IDSG2002	Island Dressing 6in x 7in Plus Adhesive Sterile	1 Each, 25 Each/Box, 4 Boxes/Case, 100 Each/Case

Non-Adherent Pads

- Absorbent pad does not adhere to the wound
- Soft, outer film and inner non-woven pad for cushion and fluid absorption
- Gentle, primary dressing ideal for lightly draining wounds
- Easy removal to minimize disruption to healing tissue
- Can be easily cut into various sizes and shapes without separating
- 1-sided and 2-sided Polypropylene options available
- Available in a range of sizes



Product Code	Product Description	Quantity
NADSG101	Non-Adherent Dressing 2in x 3in Sterile	1 Each, 100 Each/Box, 24 Boxes/Case, 2400 Each/Case
NADSG102	Non-Adherent Dressing 3in x 4in Sterile	1 Each, 100 Each/Box, 24 Boxes/Case, 2400 Each/Case
NADSG103	Non-Adherent Dressing 3in x 4in Sterile	1 Each, 50 Each/Box, 18 Boxes/Case, 900 Each/Case
NADSG104	Non-Adherent Dressing 3in x 6in Sterile	1 Each, 50 Each/Box, 15 Boxes/Case, 750 Each/Case
NADSG105	Non-Adherent Dressing 3in x 8in Sterile	1 Each, 50 Each/Box, 12 Boxes/Case, 600 Each/Case
NADSG201	Non-Adherent Dressing 3in x 8in 1-Sided Polypropylene Non Sterile	200 Each/Box, 5 Boxes/Case
NADSG202	Non-Adherent Dressing 8in x 10in 1-Sided Polypropylene Non Sterile	125 Each/Box, 4 Boxes/Case
NADSG203	Non-Adherent Dressing 3in x 8in 2-Sided Polypropylene Non Sterile	200 Each/Box, 20 Boxes/Case