

ARMOR-FLEX* STERILIZATION WRAP

with ASSURE-DRY* Technology



FLEX WITH CONFIDENCE

Built for SPD and O.R. teams managing heavier trays and tighter schedules.

ARMOR-FLEX* Sterilization Wrap with ASSURE-DRY* Technology combines exceptional strength, engineered drapability, and advanced moisture management in one high-performance sterilization wrap. Its lighter-weight, conformable design makes wrapping easier, while the innovative inner layer helps manage moisture to support consistent drying and greater confidence throughout the sterilization process.

THE GOAL: fewer reprocesses, smoother setups, and greater confidence from the first fold to the final open.



KEY PRODUCT FEATURES



Superior Strength

ARMOR-FLEX* offers balanced protection with superior breach resistance. Unlike other heavy-weight wrap that relies on a higher basis weight to achieve durability, our advanced hybrid technology delivers superior strength at a lighter weight. The result is a wrap that's not only stronger, but is also less bulky and easier to handle.



Engineered Drapability

Stop wrestling with your heavy-weight wrap. Precision engineered for conformability, **ARMOR-FLEX*** drapes naturally to simplify the wrapping process and support a smooth aseptic presentation. Greater flexibility can help reduce user fatigue, improve handling, and minimize disruptions in the SPD.




Drying Assurance

Bulky wrap can create drying challenges that could lead to wet packs, reprocessing, and costly delays. **ARMOR-FLEX*** with **ASSURE-DRY* Technology** manages moisture more effectively, providing more consistent drying and greater confidence throughout the sterilization process.

A LEGACY OF INNOVATION. A NEW STANDARD OF PERFORMANCE.


HALYARD* has been at the forefront of innovation in sterilization for over 40 years. **ARMOR-FLEX***, our latest innovation, is precision engineered with new manufacturing technology to deliver more confidence in the O.R.



WHAT IS ASSURE-DRY* TECHNOLOGY?

ASSURE-DRY* Technology sets a new standard in advanced moisture management. Its unique amber inner layer features heat-activated technology that activates during the steam sterilization cycle to optimize moisture management.

As moisture is immediately absorbed and dispersed, it evaporates more efficiently while minimizing pooling. Because the technology works independently of airflow, it promotes more consistent drying and greater confidence in the sterilization process.



ENGINEERED FOR YOUR HEAVIEST SETS

Product Code	Size	Applications / Case
ARMOR24	24in. X 24in. (60 cm X 60 cm)	96
ARMOR36	36in. X 36in. (91 cm X 91 cm)	60
ARMOR40	40in. X 40in. (101 cm X 101 cm)	40
ARMOR45	45in. X 45in. (114 cm X 114 cm)	40
ARMOR48	48in. X 48in. (121 cm X 121 cm)	40
ARMOR54	54in. X 54in. (137 cm X 137 cm)	40



HALYARD* is an Owens & Minor brand

*Registered Trademark or Trademark of O&M Halyard or its affiliates. ©2026. All rights reserved. COPY-06211



www.halyardhealth.com