

YOU'RE FIGHTING CANCER.

Wear the Battle Gear That Protects Against 52 Chemotherapy Drugs!



It's 52 vs. You.

PURPLE NITRILE-XTRA* Gloves and HALYARD* Procedure Gowns for Chemotherapy Use have been tested for resistance to 52 chemotherapy drugs.†

Arsenic Trioxide (1 mg/ml)	Cyclophosphamide (20 mg/ml)	Fludarabine (25 mg/ml)	Mitomycin (0.5 mg/ml)	Temsirolimus (25 mg/ml)
Azacitidine (Vidaza) (25 mg/ml)	Cytarabine HCl (100 mg/ml)	Fluorouracil (50 mg/ml)	Mitoxantrone (2 mg/ml)	Trastuzumab (21 mg/ml)
Bendamustine (5 mg/ml)	Cytovene (10 mg/ml)	Fulvestrant (50 mg/ml)	Oxaliplatin (2 mg/ml)	ThioTEPA (10 mg/ml)
Bleomycin Sulfate (15 mg/ml)	Dacarbazine (10 mg/ml)	Gemcitabine (38 mg/ml)	Paclitaxel (6 mg/ml)	Topotecan HCl (1 mg/ml)
Bortezomib (Velcade) (1 mg/ml)	Daunorubicin HCl (5 mg/ml)	Idarubicin (1 mg/ml)	Paraplatin (10 mg/ml)	Triclosan (1 mg/ml)
Busulfan (6 mg/ml)	Decitabine (5 mg/ml)	Ifosfamide (50 mg/ml)	Pemetrexed (25 mg/ml)	Trisenox (0.1 mg/ml)
Carboplatin (10 mg/ml)	Docetaxel (10 mg/ml)	Irinotecan (20 mg/ml)	Pertuzumab (30 mg/ml)	Vincristine Sulfate (1 mg/ml)
Carfilzomib (2 mg/ml)	Doxorubicin HCl (2 mg/ml)	Mechlorethamine HCl (1 mg/ml)	Raltitrexed (0.5 mg/ml)	Vinblastine (1 mg/ml)
Carmustine (3.3 mg/ml)†	Elence (2 mg/ml)	Melphalan (5 mg/ml)	Retrovir (10 mg/ml)	Vinorelbine (10 mg/ml)
Cetuximab (Erbix) (2 mg/ml)	Eribulin Mesylate (0.5 mg/ml)	Methotrexate (25 mg/ml)	Rituximab (10 mg/ml)	Zoledronic Acid (0.8 mg/ml)
Cisplatin (1 mg/ml)	Etoposide (20 mg/ml)			

† Testing measured no breakthrough at the Standard Breakthrough Rate of 0.1ug/cm²/minute, up to 240 minutes for gloves and 480 minutes for gowns, except for carmustine. Carmustine breakthrough occurred at 80.4 minutes on gloves and 33.3 minutes on gowns.

There is no safe level of exposure¹ to chemotherapy drugs for healthcare workers.

Is your team at risk? New, stricter guidelines — **USP <800>** — apply to *all healthcare personnel* who handle hazardous drugs.²

Don't wait: get full protection now with **HALYARD* Chemo Gear!**

PURPLE NITRILE-XTRA* Exam Gloves

Code	Description	Case Count
55090	X-Small, non-sterile, ambidextrous	500
50601	Small, non-sterile, ambidextrous	500
50602	Medium, non-sterile, ambidextrous	500
50603	Large, non-sterile, ambidextrous	500
50604	X-Large, non-sterile, ambidextrous	500

More available at halyardhealth.com.

NEW PINK UNDERGUARD* Nitrile Exam Gloves

Code	Description	Case Count
47452	X-Small, non-sterile, ambidextrous	1000
47453	Small, non-sterile, ambidextrous	1000
47454	Medium, non-sterile, ambidextrous	1000
47455	Large, non-sterile, ambidextrous	1000
47456	X-Large, non-sterile, ambidextrous	1000

More available at halyardhealth.com.

HALYARD* FACE SHIELDS

Code	Description	Case Count
41204	Full length clear visor, 9.5"	40
41205	Three-quarter length clear visor, 7.5"	40

NEW CHEMO360*

Procedure Gowns Tested for Chemotherapy Use

Code	Description	Case Count
47346	Universal, Closed back, Tape closure, Poly-coated	100
47347	XXL, Closed back, Tape closure, Poly-coated	100

HALYARD* PROCEDURE GOWNS Tested for Chemotherapy Use

Code	Description	Case Count
69606	Universal, Open back, Over-the-head, Poly-coated	100
37284	XX Large, Over-the-head, Poly-coated	100

FLUIDSHIELD* N95 RESPIRATOR MASKS

Code	Description	Case Count
46828	Small, individually bagged	210
46728	Regular, individually bagged	210
46827	Small	210
46727	Regular	210
46867	Small, with SAFETY SEAL*	210
46767	Regular, with SAFETY SEAL*	210

More available at halyardhealth.com.



Get the facts. Visit
www.GearUpForTheFight.com

Call 1-844-HALYARD (1-844-425-9273)
in the United States and Canada.

1. Centers for Disease Control and Prevention, Hazardous Drug Exposures in Health Care, www.cdc.gov/-niosh/topics/hazdrug/default.html
2. USP General Chapter <800> Hazardous Drugs – Handling in Healthcare Settings, United States Pharmacopoeial Convention, USP 40–NF 35, Second Supplement (2017), 1. www.usp.org

* Registered Trademark or Trademark of O&M Halyard or its affiliates. © 2018. All rights reserved.
* PURPLE, PURPLE NITRILE and PURPLE NITRILE X-TRA are Registered Trademarks or Trademarks of O&M Halyard or its affiliates.
THE COLOR PURPLE is a Registered Trademark or Trademark O&M Halyard or its affiliates.