

INTRODUCING HALYARD* ALGINATE WOUND DRESSINGS



PERFORMANCE. COMFORT. SAVINGS. TREAT WITH CONFIDENCE.

HALYARD* Alginate Wound Dressing is a sterile, non-woven calcium-sodium alginate dressing consisting of high strength gel alginate.

- For use on moderate-to-heavily exuding wounds and partial and full thickness wounds
- Highly absorbent, conformable and maintains structural integrity when wet
- Promote a moist wound-healing environment
- Can be used under compression

OUR NEW WOUND CARE LINE IS PART OF THE HALYARD* PORTFOLIO



INFECTION CONTROL



INCONTINENCE



SURGICAL SOLUTIONS



ADVANCED WOUND CARE

HALYARD*

Alginate Wound Dressings

EXUDATE LEVEL	PRODUCT CODE	DESCRIPTION	CASE QTY.
	50096	HALYARD* Alginate, 2 in. x 2 in. (5cm x 5cm)	100
	50097	HALYARD* Alginate, 4 in. x 4 in. (10cm x 10cm)	100
	50098	HALYARD* Alginate, 4 in. x 8 in. (10cm x 20cm)	50
	50099	HALYARD* Alginate, 1 in. x 12 in. (2.5cm x 30cm)	50
	50100	HALYARD* Alginate, 4.25 in. x 4.25 in. (10.8cm x 10.8cm)	100

EXUDATE LEVEL KEY:

 = Light

 = Moderate

 = Heavy

Partner with O&M Halyard to stay in the forefront of patient care.
For more information, please visit: www.halyardhealth.com or contact your
O&M Halyard representative **1-844-HALYARD** (1-844-425-9273).

www.halyardhealth.com  