

INTRODUCING HALYARD* XEROFORM WOUND DRESSINGS

**PERFORMANCE.
COMFORT. SAVINGS.
TREAT WITH CONFIDENCE.**

HALYARD* Xeroform Wound Dressing is a sterile, non-adherent protective dressing consisting of fine mesh absorbent gauze impregnated with USP White Petrolatum and 3% Bismuth Tribromophenate.

- Non-adherent and protective
- Covers and protects wounds from environmental contamination
- Mildly deodorizing
- Promotes a moist wound healing environment



EXUDATE LEVEL	PRODUCT CODE	DESCRIPTION	CASE QTY.
	50197	HALYARD* Xeroform Dressing, 5 in. x 9 in. (12.7cm x 22.8cm)	200
	50198	HALYARD* Xeroform Dressing, 2 in. x 2 in. (5cm x 5cm)	150
	50199	HALYARD* Xeroform Dressing, 4 in. x 4 in. (10cm x 10cm)	150
	50200	HALYARD* Xeroform Dressing, 1 in. x 8 in. (2.5cm x 20cm)	200

EXUDATE LEVEL KEY: = Light = Moderate = Heavy

Partner with O&M Halyard to stay in the forefront of patient care.
For more information, please visit: www.halyardhealth.com or contact your O&M Halyard representative **1-844-HALYARD** (1-844-425-9273).