








QUICK REFERENCE GUIDE

HALYARD* Nitrile Exam Gloves help safeguard your staff and patients from exposure to chemicals and other infectious materials. We offer a range of protection and value to deliver the best sterile or non-sterile glove for the task.

HIGH USE NITRILE GLOVES

Product / Color	Size	Code	Cuff Length	Fingertip Thickness	Palm Thickness	Packaging	USP <800> Compliant ¹	Chemo Drug Tested ²	Chemical Tested ³	Sterile Available
 PURPLE NITRILE*	XS	55080	9.5"	0.15 mm 5.9 mil	0.12 mm 4.7 mil	100 ea/bx, 10 bx/cs (XL: 90 ea/bx, 10 bx/cs)	✓	✓	✓	✓
	S	55081								
	M	55082								
	L	55083								
	XL	55084								
 STERLING*	XS	50705	9.5"	0.12 mm 4.7 mil	0.08 mm 3.1 mil	200 ea/bx, 10 bx/cs (XL: 170 ea/bx, 10 bx/cs)	✓	✓	✓	✓
	S	50706								
	M	50707								
	L	50708								
	XL	50709								
 STERLING SG*	XS	41657	9.5"	0.09 mm 3.5 mil	0.08 mm 3.1 mil	250 ea/bx, 10 bx/cs (XL: 230 ea/bx, 10 bx/cs) Available in ECO PULL* Dispensing System (47339 – 47343)	✓	✓	✓	
	S	41658								
	M	41659								
	L	41660								
	XL	41662								
 LAVENDER*	XS	52816	9.5"	0.08 mm 3.1 mil	0.07 mm 2.8 mil	250 ea/bx, 10 bx/cs (XL: 230 ea/bx, 10 bx/cs)	✓	✓	✓	
	S	52817								
	M	52818								
	L	52819								
	XL	52820								
 AQUASOFT*	XS	43932	9.5"	0.08 mm 3.1 mil	0.07 mm 2.8 mil	300 ea/bx, 10 bx/cs (XL: 250 ea/bx, 10 bx/cs)	✓	✓	✓	
	S	43933								
	M	43934								
	L	43935								
	XL	43936								

HIGH USE VINYL GLOVES

Product / Color	Size	Code	Cuff Length	Fingertip Thickness	Palm Thickness	Packaging	USP <800> Compliant ¹	Chemo Drug Tested ²	Chemical Tested ³	Sterile Available
 SYNTHETIC PLUS	S	50031	9.5"	0.09 mm 3.5 mil	0.11 mm 4.3 mil	100 ea/bx, 10 bx/cs (XL: 90 ea/bx, 10 bx/cs)				
	M	50032								
	L	50033								
	XL	50034								
 SYNTHETIC	XS	55030	9.5"	0.09 mm 3.5 mil	0.11 mm 4.3 mil	100 ea/bx, 10 bx/cs (XL: 90 ea/bx, 10 bx/cs)				
	S	55031								
	M	55032								
	L	55033								
	XL	55034								

¹ Meets the USP <800> guidelines for safe handling of hazardous drugs

² Tested in accordance with the ASTM D6978 standard

³ Tested in accordance with the ASTM F739 standard

All physical properties and product specs are target measurements; target MIL measurements may differ due to rounding.

Specialty

QUICK REFERENCE GUIDE

From proprietary nitrile formulations and specialized glove designs, to enhanced tactile sensitivity and economical options, there's a HALYARD* Medical Exam Glove that's right for every clinical task.

SPECIALTY GLOVES

Product / Color	Size	Code	Specialty Area	Cuff Length	Fingertip Thickness	Palm Thickness	Packaging	USP <800> Compliant ¹	Chemo Drug Tested ²	Chemical Tested ³	Sterile Available
PURPLE NITRILE* MAX	S	44992	Decontamination	16"	0.28 mm 11.0 mil	0.24 mm 9.4 mil	50 ea/bx 8 bx/cs	✓	✓ 128.5 > 240	C T	✓
	M	44993									
	L	44994									
	XL	44995									
PURPLE NITRILE-XTRA*	XS	55090	Chemo drug handling (overglove)	12"	0.15 mm 5.9 mil	0.12 mm 4.7 mil	50 ea/bx 10 bx/cs	✓	✓ 80.4 > 240	C T	✓
	S	50601									
	M	50602									
	L	50603									
PINK UNDERGUARD*	XS	47452	Chemo drug handling (underglove)	12"	0.12 mm 4.7 mil	0.09 mm 3.5 mil	100 ea/bx 10 bx/cs	✓	✓ 54.7 128.5	C T	✓
	S	47453									
	M	47454									
	L	47455									
STERLING NITRILE-XTRA*	XS	53137	Chemo drug handling	12"	0.12 mm 4.7 mil	0.09 mm 3.5 mil	100 ea/bx 10 bx/cs	✓	✓	✓	✓
	S	53138									
	M	53139									
	L	53140									
STERLING ZERO*	XS	47459	Manufactured without added accelerators	9.5"	0.10 mm 3.9 mil	0.07 mm 2.8 mil	200 ea/bx 10 bx/cs (XL: 170 ea/bx, 10 bx/cs)	✓	✓	✓	
	S	47460									
	M	47461									
	L	47462									
BLACK-FIRE*	XS	44755	EMS, First Responders, Law Enforcement	9.5"	0.14 mm 5.5 mil	0.10 mm 3.9 mil	150 ea/bx, 10 bx/cs (XL: 140 ea/bx, 10 bx/cs)	✓			✓
	S	44756									
	M	44757									
	L	44758									
FLEXAPRENE* GREEN	XS	44792	High tactile sensitivity procedures	9.5"	0.09 mm 3.5 mil	0.07 mm 2.8 mil	200 ea/bx, 10 bx/cs (XL: 180 ea/bx, 10 bx/cs)	✓			
	S	44793									
	M	44794									
	L	44795									
XL	44796										