

HALYARD* FLUIDSHIELD*

Surgical N95 Mask



MEDICAL GRADE. MADE IN THE AMERICAS.

Not all N95s are created equal. Make sure yours are up to healthcare standards with HALYARD* FLUIDSHIELD* Surgical N95 masks.

100% MEDICAL GRADE FROM THE LEADER IN TOTAL PPE FOR HEALTHCARE

While some N95s approved for emergency use are industrial grade, 100% of HALYARD* N95s are **FDA-cleared medical devices** for use in healthcare settings.

MADE IN THE AMERICAS, FROM FABRIC TO FINISH

Unlike suppliers who rely on sources from Asia, all FLUIDSHIELD* Surgical N95s are **manufactured here** in the Americas.

HIGH-LEVEL SPLASH PROTECTION

Our FLUIDSHIELD* Level 3 respirators provide your staff with the highest-recognized level of protection from splash.¹

BREATHABILITY AND FIT DESIGNED SPECIFICALLY FOR MEDICAL STAFF

The FLUIDSHIELD* duckbill breathing chamber is more than **twice as large** as the leading surgical N95, and exceeds NIOSH standards for breathability.

Because proper snug fit is critical², HALYARD* N95s come in a range of sizes, with malleable nose wire to adjust fit. They have strong elastic straps that are securely **bonded—not stapled**—to the mask.

RELIABLE SUPPLY EVERY DAY

With manufacturing investment and end-to-end control of our **entire supply chain**, we have been able to ramp up N95 production to meet customer allocation levels even during the pandemic.




¹ ASTM F2100 Level 3 splash protection tested per ASTM F1862 at 160 mmHg pressure.

² Each wearer must be fit-tested for the brand and model of N95 they will be using

BE READY FOR THE FUTURE WITH HALYARD* MEDICAL GRADE N95 RESPIRATORS

It makes sense to stockpile the PPE that you will need to help protect your staff, whatever the future brings. HALYARD* FLUIDSHIELD* Surgical N95s are:

- NIOSH approved to filter out at least 95% of very small (<0.3micron) particles, including bacteria and viruses.
- FDA-cleared Class 2 medical devices
- SO SOFT* lined for comfort
- Folds flat for compact storage
- Meets — and exceeds — NIOSH breathability standards



**MADE IN THE AMERICAS
FOR SUPPLY RELIABILITY**

FLUIDSHIELD* Surgical N95

CODE	DESCRIPTION	SIZE	BOXES/CASE	CASE COUNT
46727	Level 3 splash-resistant, NIOSH & FDA 510k Approved, latex-free	REGULAR	6	210
46827	Level 3 splash-resistant, NIOSH & FDA 510k Approved, latex-free	SMALL	6	210



Ask about our HALYARD* Qualitative Fit Test Kit

Contains all the materials and instructions you need to conduct Fit Testing of N95 Respirators. Each kit contains supplies for 75-100 fit tests. Order code 47950.



Knowledge Network* Accredited Education
In-Service Training and Technical Support
Credentialed Sales Representatives
Tools & Best Practices
Clinical Research
Commitment to Excellence

Talk to your Halyard representative
Or visit www.halyardhealth.com