

MEDICAL ACTION MINOR PROCEDURE KITS

Did you know NDC is now stocking all of these Medical Action Minor Procedure Kits? Check the NDC website for pricing.

IV START

#71952 – 100/CS	Qty/Kit	UOM
Dressing – 2-3/8" X 2-3/4"	1	EA
Gauze – 2" X 2" 4-Ply NW	2	EA
Chloraprep – 1 mL	1	EA
Tape – Transparent 0.75" X 18"	1	EA
Tourniquet – 0.75" X 18" Rolled Blue	1	EA

CENTRAL LINE DRESSING CHANGE

#262834 – 20/CS	Qty/Kit	UOM
Tegaderm – 4" X 4.75" (3M 1626WNS)	1	EA
Label – Dressing Change 1.6" X 1"	1	EA
Gauze – 4" X 4" 4-Ply NW	2	EA
Tape – Transparent Strip 1" X 4"	1	EA
Alcohol Swab – 3/pk	1	EA
Chloraprep – 3 mL	1	EA
Glove – Nitrile Med Blue Walleted	1	PR
Wrap – 17" X 19" Blue/White	1	EA
Mask – Ear loop Folded	1	EA

CENTRAL LINE DRESSING CHANGE

#262808 – 20/CS	Qty/Kit	UOM
Tray – Clear	1	EA
Tegaderm – 4" X 4.75" (3M 1626WNS)	1	EA
Label – Dressing Change 1.6" X 1"	1	EA
Gauze – 4" X 4" 8-Ply	2	EA
Gauze – 2" X 2" 4-Ply NW	2	EA
Chloraprep – 3 mL	1	EA
Tape – Transparent 0.75" X 18"	1	EA
Alcohol Swab 1	1	EA
Applicator – Cotton Tip 6"	1	EA
Glove – Nitrile Med Blue Walleted	1	PR
Wrap – 17" X 19" Blue/White	1	EA
Mask – Ear loop Folded	1	EA

STAPLE REMOVAL

#68257 – 50/CS	Qty/Kit	UOM
Staple Remover – 4"	1	EA
Gauze – 3" X 3" 12-Ply	1	EA

SUTURE REMOVAL – LITTAUER SCISSOR

#M2634 – 50/CS	Qty/Kit	UOM
Scissor – 4.5" Iris (with point protector)	1	EA
Forcep – 4.75" Adson Serrated (with point protector)	1	EA
Gauze – 2" X 2" 8-Ply	2	EA
PVP Iodophor Prep Pad Med	1	EA
Alcohol Med Prep Pad Med	1	EA

SUTURE REMOVAL – LITTAUER SCISSOR (no antiseptic)

#4131 – 50/CS	Qty/Kit	UOM
Scissor – 4.5" Wire Littauer (with point protector)	1	EA
Forcep – 5" Thumb Wire Serrated	1	EA
Gauze – 2" X 2" 8-Ply	1	EA

LACERATION TRAY (no needles)

#61280 – 20/CS	Qty/Kit	UOM
Med Cup – 60 cc Gamma Safe	2	EA
Inst: Bundle – Forcep/Scissor/Hemo/NH	1	EA
Gauze – 4" X 4" 8-Ply	10	EA
Gauze – 2" X 2" 8-Ply	5	EA
Towel – Blue Blotting 4-Ply	2	EA
Drape – 18" X 26" with Fen 3" - White Tissue/Poly	1	EA
Towel – 13" X 18" 2-Ply Folded	1	EA

LACERATION TRAY (with needles)

#61281 – 20/CS	Qty/Kit	UOM
Med Cup – 60 cc Gamma Safe	2	EA
Inst: Bundle – Forcep/Scissor/Hemo/NH	1	EA
Gauze – 4" X 4" 8-Ply	10	EA
Gauze – 2" X 2" 8-Ply	5	EA
Syringe – 12 cc L/L	1	EA
Needle – 25 Ga X 0.625" Orange	1	EA
Needle – 18 Ga X 1.5" Pink	1	EA
Needle – 22 Ga X 1.5" Dark Grey	1	EA
Needle – 27 Ga X 1.25" Grey	1	EA
Towel – Blue Blotting 4-Ply	2	EA
Drape – 18" X 26" with Fen 3" - White Tissue/Poly	1	EA
Towel – 13" X 18" 2-Ply Folded	1	EA

The complete line of Medical Action Minor Procedure Kits & Trays are available through NDC.