



PRIMO Defender

Defender offloading heel boot is an economical and simple way to offload the heel. Filled with top quality soft hollow core fiber, the heel comfortably floats above the bed protecting against heel pressure injuries.



Maintain correct foot placement with interior soft micro-suede fabric and high sides to wrap leg

Support the prevention of skin shearing with durable outer fabric that repels fluid and slides easily



Heel is completely offloaded with a large heel hole



*Standard Size, Purple Floats heel by using top quality soft hollow core fiber



Secure a snug fit with two strap system

PRIMO Defender			
Item #	Size and Calf Measurement	Color	Case QTY
HB100S	Small - Less than 12.5"	Black	12 Each
HB100	Standard - 12.5" - 17"	Purple	12 Each
HB100B	Bariatric - Greater than 17"	Navy	10 Each

PRIMO Guardian Heel Protector			
Item #	Size and Calf Measurement	Color	Case QTY
HP201N	Standard - 12.5" - 17"	NAVY/GRAY	12 Each



Bariatric



Small