



PRIMO HeelCheck

HeelCheck offloading heel boot is designed for an easy heel check. Clinically designed to increase compliance, plantar panels open to assess the skin, check pulses, and cool the foot.

Prevent foot from dropping with adjustable "Y" straps that keep the foot at a 90° angle



Maintain correct foot placement with interior soft micro-suede fabric and high sides to wrap leg

Tapered foam wedge combined with soft pillow fiber constantly floats the heel and mimics the curve of the calf to reduce hyper extension of the knee



SCD compatible with large openings on both sides to route SCD tubes away from skin

Simplify boot application by using a numbered strap system

FDA REQUIRES HEEL OFFLOADING BOOTS ARE SINGLE PATIENT USE ONLY
*It is recommended to remove entire boot for proper skin assessment.

PRIMO HeelCheck

Support the prevention of skin shearing with durable outer fabric that repels fluid and slides easily

Plantar panels open for an easy heel check



Maintain correct foot placement with interior soft micro-suede fabric and high sides to wrap leg



Heel is completely offloaded with a large heel hole

PRIMO HeelCheck			
Item #	Size and Calf Measurement	Color	Case Qty
HCS	Small - Less than 12.5"	Green/Gray	13 Each
HC1	Standard - 12.5" - 17"	Purple/Navy	13 Each
HCB	Bariatric - Greater than 17"	Navy/Gray	10 Each
ARW01	Anti-Rotational Wedge fits all HeelCheck Sizes	Black	24 Each

PRIMO HeelCheck with Anti-Rotational Wedge		
Item #	Size and Calf Measurement	Case Qty
HCSW	Small - Less than 12.5"	8 Each
HC1W	Standard - 12.5" - 17"	8 Each
HCBW	Bariatric - Greater than 17"	8 Each



Aids in prevention of internal/external rotation of the leg by attaching an anti-rotational wedge to either side of the boot (Optional Feature)