



PRIMOBoot™

PRIMOBoot, clinically designed with an internal tapered foam wedge, consistently offloads the heel providing the ultimate protection against heel pressure injuries.

Heel Lift System

The **only** heel offloading boot on the market that has a **tapered foam wedge combined with soft pillow fiber**. The PRIMOBoot heel lift system **ensures continuous flotation of the heel** for treatment and prevention of a pressure point which can lead to a pressure injury.



Prevent foot from dropping with adjustable "Y" straps that keep the foot at a 90° angle

Maintain correct foot placement with interior soft micro-suede fabric and high sides to wrap leg

Patented tapered foam wedge combined with soft pillow fiber constantly floats the heel and mimics the curve of the calf to reduce hyper extension of the knee

Support the prevention of skin shearing with durable outer fabric that repels fluid and slides easily

Simplify boot application by using a numbered strap system

SCD compatible with large openings on both sides to route SCD tubes away from skin



FDA REQUIRES HEEL OFFLOADING BOOTS ARE SINGLE PATIENT USE ONLY



For more information, contact your Owens & Minor sales representative.

©2022, 2023 All Rights Reserved. ASP Global 008565



PRIMOBoot™

PRIMOBoot				PRIMOBoot with Anti-Rotational Wedge	
Item #	Case QTY	Description	Size and Calf Measurement	Item #	Case QTY
BABBS	13 Each	Black/Gray	Small - 9" - 12.5"	BABBSW	8 Each
BAB1000	13 Each	Black/Khaki	Standard - 12.5" - 17"		
BAB11000	13 Each	Purple/Navy	Standard - 12.5" - 17"	BAB11000W	8 Each
BABB	10 Each	Navy/Gray	Bariatric - Greater than 17"	BABBW	8 Each
BABSS01	12 Each	Diabetic Boot - Lined with synthetic sheepskin. Excellent for vascular patients.	Standard - 12.5" - 17"	BABSS01-W	8 Each
ARW01	24 Each	Anti-Rotational Wedge	Fits all PRIMOBoot™ sizes except Small and Child		



Heel is completely offloaded with a large heel hole



Aids in prevention of internal/external rotation of the leg by attaching an anti-rotational wedge to either side of the boot (Optional Feature)



PRIMOBoot Child is specifically designed with a smaller fit for children.

Burn Boot & Fall Risk Boots Available

Patent No:
US 8,834,396 B2