




# HALYARD\* PUREZERO\* HG3 White Sterile Nitrile Gloves For the Critical Cleanroom Environment For Industrial Use Only

## TECHNICAL DATA SHEET

### Description

HALYARD\* **PUREZERO\*** HG3 White Sterile Nitrile Cleanroom Gloves are designed for critical cleanroom environments such as pharmaceutical and biotechnology cleanroom manufacturing as well as sterile compounding cleanroom applications. These hand specific gloves are clean processed (washed repeatedly in deionized water) to ensure consistent control of low particles, extractables and endotoxin levels, and are recommended for use in ISO Class 3 or higher and Grade A/B cleanrooms. Because HALYARD\* **PUREZERO\*** HG3 Sterile Cleanroom Gloves are made with an **accelerator-free**<sup>1</sup> nitrile polymer, there is a reduced risk of allergies and skin irritation associated with accelerator chemicals in other nitrile gloves.



Manufactured at OUR OWN  
Safeskin Facility in Thailand.

### Cleanliness Properties

Max Particle Count (>0.5 µm)	<950 Particles/cm <sup>2</sup>	IEST RP-CC005
<b>Max Endotoxin Level</b>	<20 EU	
<b>Ionic Content (Extractable ions)</b>	<b>Max Level (ug/g)</b>	IEST RP-CC005
Calcium	50	
Chloride	35	
Magnesium	5	
Nitrate	20	
Potassium	5	
Sodium	10	
Sulfate	10	
Zinc	7	
Ammonium	5	

### Physical Properties

AQL	1.0
Sterile	✓
Hand Specific Pairs	✓
Smooth Grip	✓
Textured Fingertips and Palms	✓
Accelerator-Free <sup>1</sup>	✓
Latex-Free	✓
Powder-Free	✓
Silicone-Free	✓
Static Dissipative in Use	✓
Tensile Strength <sup>2</sup>	35 MPa (Target)
Ultimate Elongation <sup>2</sup>	600%
Sterility Assurance Level (SAL)	10 <sup>-6</sup>
Shelf Life	3 Years

# TECHNICAL DATA SHEET

## Glove Dimensions

	6.0	6.5	7.0	7.5	8.0	8.5	9.0	10.0
<b>Glove Length (inch/cm)</b>	12"/30.5	12"/30.5	12"/30.5	12"/30.5	12"/30.5	12"/30.5	12"/30.5	12"/30.5
<b>Width of Palm (mm)</b>	80	87	94	98	109	114	120	128
<b>Middle finger length (mm)</b>	70.8	72.4	80.6	83.9	86.0	90.1	92.7	96.2
<b>Finger Tip Thickness</b>	0.16 mm (6.3 mil)	0.16 mm (6.3 mil)	0.16 mm (6.3 mil)	0.16 mm (6.3 mil)	0.16 mm (6.3 mil)	0.16 mm (6.3 mil)	0.16 mm (6.3 mil)	0.16 mm (6.3 mil)
<b>Palm Thickness</b>	0.13 mm (5.1 mil)	0.13 mm (5.1 mil)	0.13 mm (5.1 mil)	0.13 mm (5.1 mil)	0.13 mm (5.1 mil)	0.13 mm (5.1 mil)	0.13 mm (5.1 mil)	0.13 mm (5.1 mil)
<b>Cuff Thickness</b>	0.10 mm (3.9 mil)	0.10 mm (3.9 mil)	0.10 mm (3.9 mil)	0.10 mm (3.9 mil)	0.10 mm (3.9 mil)	0.10 mm (3.9 mil)	0.10 mm (3.9 mil)	0.10 mm (3.9 mil)

## Packaging Data

Triple layer packaging (poly pouch and poly bag plus case liner)

200 hand specific pairs per case: 2 gloves/poly wallet & pouch x 20 sealed pouches per bag x 10 PE bags per lined carton

Packaged in ISO Class 5 Cleanroom

## Quality & Regulatory Standards

Compliant to these regulatory standards:

- ISO 9001
- ISO 13485
- ISO 14001

Compliant to these food handling regulatory standards:

- FDA 21 CFR 177-2600
- FDA 21 CFR 180.22
- Commission Regulation (EU) No 10/2011
- CE 2797 PPE Category III according to Regulation (EU) 2016/425 EEC
- EN 16523-1:2015+A1:2018
- EN ISO 374-2:2019
- EN ISO 374-4:2019
- EN ISO 374-5:2016 Bacteria, Fungi, and Virus Protection
- EN ISO 374-1:2016+A1:2018/ Type B
- EN ISO 21420:2020 Dexterity Classification: Level 5


Compliant with the REACH regulation

RoHS 2011/65/EU Annex II

Static Dissipative in Use

Tested per ANSI/ ESD SP 15.1

Meets USP <800> Guidelines



## Additional Glove Information

Recommended for use in ISO Class 3 or higher and Grade A/B cleanrooms.

Tested against 29 chemicals, 14 chemotherapy drugs and Fentanyl.

Made in Thailand

Declaration of Conformity (DoC), Certificates of Analysis (COA) and Certificates of Processing (COP) for every production lot available online at [halyardhealth.com/information](http://halyardhealth.com/information)

## Storage Instructions

HALYARD\* PUREZERO\* Nitrile Gloves shall be stored in conditions where the product is kept dry (away from moisture), away from direct sunlight, away from sources of heat, and away from radioactive sources.

## Ordering Information

**HALYARD\* PUREZERO\* HG3 WHITE STERILE NITRILE GLOVES, HAND SPECIFIC, SMOOTH GRIP, TEXTURED FINGERTIPS AND PALMS**

Size	Code	Size	Code
6.0	CLN323260	8.0	CLN323280
6.5	CLN323265	8.5	CLN323285
7.0	CLN323270	9.0	CLN323290
7.5	CLN323275	10.0	CLN323210

For additional information or samples, contact your local distributor or visit [www.purezerogloves.com](http://www.purezerogloves.com)

1 Not formulated with these commonly used vulcanizing chemicals: Sulfur, Thiurams, Thioxoles, Guanidines and Carbamates.

2 Tested per ASTM D412

This fact sheet has been created using the most recent information. In the interest of continuous improvement, the characteristics of the products may change without prior notice.

\*Registered Trademark or Trademark of O&M Halyard or its affiliates. ©2020, 2021, 2023. All rights reserved. 008640 COPY-04709 September 2023