



# HALYARD

*Part of the Owens & Minor Family*

# DELIVERING CONFIDENCE.

## PHYSICIAN OFFICE CATALOG

- GLOVES • FACIAL PROTECTION • GOWNS • DRAPES
- STERILIZATION WRAP • MINOR PROCEDURE KITS



# You have to be able to trust your protection. With Halyard, you can.

## Halyard products are clinically superior and tested to the highest standards.

There is too much at stake - for you, your staff and your patients.

With outstanding protection, proven performance and world-class expertise, you can have total confidence in every glove, mask, gown, drape, wrap and procedure kit.

- **Customer Service** from a team that is experienced, hospital-trained and credentialed
- **Product Availability** due to vertically integrated manufacturing and distribution
- **Healthcare Expertise** that allows us to make knowledgeable recommendations

## MEDICAL EXAM GLOVES

HALYARD\* Nitrile Exam Gloves help safeguard you, your staff and patients from exposure to bloodborne pathogens, chemicals and other infectious materials. We offer a range of protection and value to deliver the best sterile or non-sterile glove for the task.

GLOVE	9.5"	12"	16"	FLUID EXPOSURE	AMBIDEXTROUS + BEADED CUFF	STERILE OPTION AVAILABLE	CHEMO DRUGS TESTED	COMMON CHEMICALS TESTED	FENTANYL RATED 240 MINS	ECO PULL* OPTION AVAILABLE
AQUASOFT*	✓			Low	✓		12	7		
LAVENDER*	✓			Low to Moderate	✓		12	20		
STERLING SG*	✓			Low to High	✓		12	10		✓
STERLING*	✓	✓		Low to High	✓	✓	12	40		
PURPLE NITRILE*	✓			Moderate to High	✓	✓	29	41	✓	
PURPLE NITRILE-XTRA*		✓		High	✓	✓	52	45	✓	
PURPLE NITRILE* MAX			✓	High	✓		52	25		
PINK UNDERGUARD*		✓		High	✓		52	6		

### AQUASOFT\* NITRILE EXAM GLOVES



#### Value, comfort, and protection

- Ultra-thin with the tactile sensitivity comparable to latex without the risk of latex skin irritation and allergies
- 67% less storage space required with 300 count dispenser box<sup>1</sup>

43932 XS  
43933 S  
43934 M  
43935 L  
43936 XL

XS - L 300 gloves/box | 10 boxes/case  
XL 250 gloves/box | 10 boxes/case

### LAVENDER\* NITRILE EXAM GLOVES



#### For procedures with low to moderate risk of exposure to fluid and abrasion

- Thinner and lighter than standard nitrile gloves with high tactile sensitivity
- #1 branded nitrile glove in the U.S. dental market<sup>2</sup>

52816 XS  
52817 S  
52818 M  
52819 L  
52820 XL

XS - L 250 gloves/box | 10 boxes/case  
XL 230 gloves/box | 10 boxes/case



## STERLING SG\* NITRILE EXAM GLOVES

**For procedures with low-to-high risk of exposure to fluid and abrasion**

- Cleared for use in chemotherapy
- Flexible with high tactile sensitivity
- Non-sterile
- Ambidextrous, textured fingertips
- Beaded cuffs

41657	XS
41658	S
41659	M
41660	L
41662	XL

XS - L	250 gloves/box   10 boxes/case
XL	230 gloves/box   10 boxes/case



## STERLING\* NITRILE EXAM GLOVES

**For procedures with low-to-high risk of exposure to fluid and abrasion**

- Cleared for use in chemotherapy
- Flexible with high tactile sensitivity
- Non-sterile
- Ambidextrous, textured fingertips
- Beaded cuffs

50705	XS
50706	S
50707	M
50708	L
50709	XL

XS - L	200 gloves/box   10 boxes/case
XL	170 gloves/box   10 boxes/case



## PURPLE NITRILE\* EXAM GLOVES

**Provides maximum barrier protection**

- Heavy weight, ideal for moderate to high risk settings and fluid exposure
- 9.5" exam glove

55080	XS
55081	S
55082	M
55083	L
55084	XL

XS - L	100 gloves/box   10 boxes/case
XL	90 gloves/box   10 boxes/case



## PURPLE NITRILE-XTRA\* EXAM GLOVES

**Provides maximum barrier protection**

- Heavy weight, ideal for moderate to high risk settings and fluid exposure
- 12" exam glove

55090	XS
50601	S
50602	M
50603	L
50604	XL

XS - L	50 gloves/box   10 boxes/case
--------	-------------------------------



## PURPLE NITRILE\* MAX EXAM GLOVES

**Longer 16" cuff for added protection, disposable**

**Provides maximum barrier protection**

- Our heaviest weight (20 grams) for added confidence when handling instruments and decontamination chemicals
- Thicker Fingertips
- Fully Textured forexcellent wet and dry grip
- Not made with natural rubber latex

44992	S
44993	M
44994	L
44995	XL

S - XL	50 gloves/box   8 boxes/case
--------	------------------------------



## PINK UNDERGUARD\* NITRILE EXAM GLOVES

**Coated exterior of glove reduces friction, making donning a second pair of gloves easy**

- Tested with ASTM D6978, approved for use with 52 chemotherapy drugs
- Extended 12" cuff providing additional coverage under the cuff of the chemo gown
- Pink color provides breach detection when wearing contrasting PURPLE NITRILE\* or STERLING\* gloves as outer layer

47452	XS
47453	S
47454	M
47455	L
47456	XL

100 gloves/box   10 boxes/case
--------------------------------

# FACIAL PROTECTION

Halyard is the leader in facial protection and our Face Masks are the choice for more healthcare professionals where comfort and protection are needed in a variety of procedures. We offer a full range of styles, features and levels of protection. From standard, to ASTM LEVEL 3 fluid resistant masks, particulate filter respirators, eyewear and face shields, you can choose the right mask for the task.



## DISNEY® CHILD'S FACE MASK

32856

The only FDA Approved and size-tested for ages 4-12

- SO SOFT\* Earloops
- Exclusively from Halyard
- Formable nose wire fits small faces
- Not made with natural rubber latex

75 masks/box  
10 boxes/case



## FLUIDSHIELD\* LEVEL 3, FOG-FREE PROCEDURE & SURGICAL MASKS, ASTM F2100-11

Protects from fluid splatter and splashing

- Meets ASTM F2100-11 Level 3 standards for fluid resistance, filtration and breathability
- Pleat-style available with earloops or ties



**28800** FLUIDSHIELD\* Level 3 Procedure Mask with Earloops and **Wraparound Visor**, Foam Band, Orange

**28804** FLUIDSHIELD\* Level 3 Surgical Mask with Ties and **Wraparound Visor**, Foam Band, Orange

25 masks/box  
4 boxes/case



**00146** FLUIDSHIELD\* Level 3 Procedure Mask with Earloops, **Wraparound Splashguard\* Visor**, Blue

**00148** FLUIDSHIELD\* Level 3 Procedure Mask with Earloops, Foam Band, Blue

**28797** FLUIDSHIELD\* Level 3 Procedure Mask with Earloops, Foam Band, Orange

**28820** FLUIDSHIELD\* Level 3 Surgical Mask with Ties, Foam Band, Orange

40 masks/box  
10 boxes/case

50 masks/box  
6 boxes/case



## FLUIDSHIELD\* LEVEL 3, FOG-FREE PROCEDURE & SURGICAL MASKS WITH SO SOFT\* LINING, ASTM F2100-11

Protects from fluid splattering and splashing

- Meets ASTM F2100-11 Level 3 standards for fluid resistance, filtration and breathability
- Pleat-style available with earloops or ties
- Features SO SOFT\* bicomponent fabric for soft comfort



**47107** FLUIDSHIELD\* Level 3 Procedure Mask, Foam Band, Earloops, Orange

**48207** FLUIDSHIELD\* Level 3 Surgical Mask, Foam Band, Ties, Orange

**47147** FLUIDSHIELD\* Level 3 Procedure Mask, **Wraparound Visor**, Foam Band, Earloops, Orange

**48247** FLUIDSHIELD\* Level 3 Surgical Mask, **Wraparound Visor**, Foam Band, Ties, Orange

40 masks/box | 10 boxes/case

50 masks/box | 6 boxes/case

25 masks/box  
4 boxes/case





## LEVEL 2, FOG-FREE PROCEDURE & SURGICAL MASKS, ASTM F2100-11

For protection from moderate risk of fluid splash and spray

- Meets ASTM F2100-11 Level 2 standards for fluid resistance, filtration and breathability
- Pleat-style available with earloops or ties



<b>62113</b>	FLUIDSHIELD* Level 2 Surgical Mask with Earloops, Blue	50 masks/box   6 boxes/case
<b>62114</b>	FLUIDSHIELD* Level 2 Surgical Mask with Ties, <b>Wraparound Visor</b> , Blue	25 masks/box   4 boxes/case
<b>62115</b>	FLUIDSHIELD* Level 2 Procedure Mask with Ties, Blue	50 masks/box   10 boxes/case
<b>62116</b>	FLUIDSHIELD* Level 2 Procedure Mask, <b>Wraparound Visor</b> , Earloops, Blue	25 masks/box   4 boxes/case



## LEVEL 1 PROCEDURE & SURGICAL MASKS, ASTM F2100-11

Meets ASTM F2100-11 Level 1 standards for fluid resistance, filtration and breathability

- Assured fluid protection for routine patient care when *low level of exposure to fluids expected*
- Pleat-style with inner layer of SO SOFT\* fabric for comfort



<b>25867</b>	FLUIDSHIELD* Level 1 Procedure Mask with SO SOFT* Lining, Earloops, Yellow	50 masks/box
<b>25868</b>	FLUIDSHIELD* Level 1 Procedure Mask with SO SOFT* Lining, Earloops, Lavender	10 boxes/case
<b>25869</b>	FLUIDSHIELD* Level 1 Procedure Mask with SO SOFT* Lining, Earloops, Blue	
<b>28802</b>	FLUIDSHIELD* Level 1 Surgical Mask with SO SOFT*Lining, Ties, Blue	50 masks/box
<b>28806</b>	FLUIDSHIELD* Level 1 <b>Fog-Free</b> Surgical Mask with SO SOFT*Lining, Green	6 boxes/case



## STANDARD PROCEDURE MASKS

Standard masks for routine patient exams – no fluid exposure expected

- Pleat-style with earloops



<b>49700</b>	Value Procedure Mask, Earloops, Blue	50 masks/box 10 boxes/case
<b>47080</b>	Standard Procedure Mask, Earloops, Blue	
<b>47085</b>	Standard Procedure Mask, Earloops, Green	
<b>47090</b>	Standard Procedure Mask, Earloops, White	
<b>47117</b>	Standard Procedure Mask, Earloops, Yellow	
<b>62362</b>	<b>Fog-Free</b> Procedure Mask, Foam Band, Earloops, Blue	
<b>62363</b>	So Soft* Procedure Mask, <b>Fog-Free</b> , Earloops, White	



## FACE SHIELDS

Protects from splash, splatter and spray

- Full or 3/4 length
- Fog-resistant
- Foam band
- Wear with or without glasses

<b>41204</b>	GUARDALL SHIELD*, Full Length Visor	40 shields/box
<b>41205</b>	GUARDALL SHIELD*, 3/4 Length Visor	1 box/case
<b>SV50A</b>	SAFEVIEW* Eyewear, May be worn over prescription glasses, assembled, assorted colors	5/box 10 boxes/case

# PROTECTIVE APPAREL

Halyard offers a full line of Protective Apparel for both staff and patients. As with our other infection prevention products, we offer apparel options to provide the best value for the task.



## HALYARD BASICS\* NON-REINFORCED SURGICAL GOWN

99284 | 20/case

HALYARD BASICS\* Surgical Gowns are designed to meet established industry performance standards for lower-fluid, shorter procedures.

- SMS fabric for barrier protection
- Raglan sleeves, adjustable hook-and-loop closures
- Also available in Extra-Large



## SPECTRUM\* NON-REINFORCED SURGICAL GOWN

46865 | 28/case

- Meets AAMI 3 Liquid Barrier Standard requirements for surgical gowns
- For procedures with low to moderate fluid contact
- Low lint, abrasion resistant
- Available in sizes M-3XL



## TRI-LAYER ISOLATION GOWN

69981 | 100/case

- Meets AAMI Level 2 requirements
- Medium-weight SMS fabric
- Available in Extra-Large.



## PATIENT GOWN

69766 | 100/case

- Single-use, full coverage gown
- Meets NFPA Class 1 standard for flame resistance
- Waist tie
- Universal size

# COLD THERAPY

Cold Packs – General, Specialty and Instant require minimal setup time with virtually no cleanup compared to ice-filled gloves, plastic bags or other makeshift ice bags.



## STAY-DRY\* ICE PACKS

33500

- 3-layer construction minimizes condensation
- Soft outer cover
- Clip closure



## INSTANT COLD PACKS

59688

Immediate cold therapy when ice is not available

- Simply squeeze and shake to activate
- Consistently maintains cold temperature for 20 minutes
- Convenient and portable



# MINOR PROCEDURE KITS

## SUTURE REMOVAL TRAYS

Choose from our wide selection of high quality, cost effective suture removal trays. We offer stainless steel, wire form, and plastic grade instruments.

50 trays/case

	ECONOMY	STANDARD	
TRAY NUMBER	61531	61537	68283
COMPONENT			
ALCOHOL MED PREP PAD			1
PVP LODOPHOR PREP PAD			1
SCISSORS	WF	SS	SS
FORCEP	WF	SS	SS
GAUZE 2X2		2	2
GAUZE 3X3	1		

WF- Wire Framed SS- Stainless Steel



## SINGLE STERILE INSTRUMENTS

Our single sterile instruments can be used for a variety of applications. Choose from our most popular items below or let your sales representative recommend the right solution for you.



Instruments 20/Box

PRODUCT ID	DESCRIPTION
56302	Scissor 5.5" OR Sharp/Blunt
56304	Scissor 4.5" Iris



# STERILIZATION WRAP

Halyard is the #1 provider of sterilization wrap in the U.S. - providing nearly 75% of all wrap used. Available in multiple sizes and grades from H100 to H650. All HALYARD\* Sterilization Wrap features POWERGUARD\* Technology for enhanced microbial barrier protection.



## QUICK CHECK\* STERILIZATION WRAP


**TWO DISTINCT LAYERS. ONE CLEAR CHOICE.** With two contrasting layers - translucent and opaque - bonded together, inspecting trays is easier and faster than ever! QUICK CHECK\* Sterilization wrap is available in multiple sizes and grades: H100, H200, H300, H400, H500 and H600.

## SMART-FOLD\* HIGH PERFORMANCE STERILIZATION WRAP

SMART-FOLD\* High Performance Wrap is 3 times more tear-resistant with reinforced zones around tray corners and edges. Using SMART-FOLD\* can minimize or eliminate costly accessories, streamline tray wrapping and can ultimately help reduce reprocessing costs. H650 is available in 5 sizes.

# DRAPES

Ask about our complete line of drapes that cover an extensive range of procedures.



## SURGICAL HALF DRAPE 60" X 44"

89101

- A versatile option from our extensive line of multi-procedural drapes
- Maximum rating for flame resistance
- Resistant to tearing, abrasion and fluid strikethrough



KNOWLEDGE NETWORK\* Accredited Education  
In-Service Training and Technical Support  
Credentialed Sales Representatives  
Tools & Best Practices  
Clinical Research  
Commitment to Excellence

**For more information visit**  
**[halyardhealth.com](http://halyardhealth.com)**

<sup>1</sup> vs. 3 boxes of 100 count gloves

<sup>2</sup> Strategic Data Marketing, 2019 Dental Distributor Data for LAVENDER\* branded gloves.

\*Registered Trademark or Trademark of O&M Halyard or its affiliates. ©2019. All rights reserved.

\*PURPLE and PURPLE NITRILE are Registered Trademarks or Trademarks of O&M Halyard or its affiliates.

©2019. All rights reserved.

The COLOR PURPLE, COLOR GRAY, COLOR LAVENDER and COLOR ORANGE, are Registered Trademarks

or Trademarks of O&M Halyard or its affiliates. ©2019. All rights reserved.

The COLOR PINK is a Registered Trademark owned by Colur World, LLC.

©2005 Colur World LLC. All rights reserved.

Disney Elements ©DISNEY

008047 COPY-03465

