

HALYARD* PUREZERO* NITRILE CLEANROOM GLOVES

HALYARD* offers an accelerator-free¹ portfolio of three non-sterile and two sterile cleanroom gloves to address the needs of pharmaceutical, medical device, microelectronics and semi-conductor manufacturing industries.

GLOVE SELECTION GUIDE

	Description	Designed for	Max Particle Count	ISO Class	Finish	Fingertips	Double Donning	Size Range	Fingertip Thickness	Case Count
NON-STERILE	HALYARD* PUREZERO* HG3 White Nitrile Gloves	Semiconductor, Pharmaceutical, Medical Device Manufacturing	<950	ISO Class 3 or higher and Grade A/B/C/D cleanrooms	Tacky Grip	Textured Fingertips	Recommended for Outer Glove	XS-XL	.16mm (6 mil)	1000/cs Ambidextrous
	HALYARD* PUREZERO* HG3 Light Blue Nitrile Gloves	Semiconductor, Pharmaceutical, Medical Device Manufacturing	<1200	ISO Class 3 or higher and Grade A/B/C/D cleanrooms	Tacky Grip	Textured Fingertips	Recommended for Outer Glove	XS-XL	.10mm (4 mil)	1500/cs Ambidextrous
	HALYARD* PUREZERO* HG3 White SGX* Nitrile Gloves	Semiconductor, Pharmaceutical, Medical Device Manufacturing	<950	ISO Class 3 or higher and Grade A/B/C/D cleanrooms	Smooth Grip	Textured Fingertips	Outer or Under	XS-XL	.16mm (6 mil)	1000/cs Ambidextrous
STERILE	HALYARD* PUREZERO* HG3 Light Blue Sterile Nitrile Gloves	Pharmaceutical, Biotechnology, Sterile Compounding, Aseptic Processing	<1200	ISO Class 3 or higher and Grade A/B/C/D cleanrooms	Smooth Grip	Textured Fingertips and Palms	Outer or Under	6.0 6.5 7.0 7.5 8.0 8.5 9.0 10.0	.10mm (4 mil)	300 pairs/cs Hand Specific
	HALYARD* PUREZERO* HG3 White Sterile Nitrile Gloves	Pharmaceutical, Biotechnology, Sterile Compounding, Aseptic Processing	<950	ISO Class 3 or higher and Grade A/B/C/D cleanrooms	Smooth Grip	Textured Fingertips and Palms	Outer or Under	6.0 6.5 7.0 7.5 8.0 8.5 9.0 10.0	.16mm (6 mil)	200 pairs/cs Hand Specific
Applies to all PUREZERO* Gloves:		Accelerator-Free ¹ Static Dissipative in Use ² Low Dermatitis Potential Compliant with the REACH regulation		AQL 1.0	CE 2797 (PPE Cat. III)	ISO 374-5:2016 VIRUS	ISO 374-1:2016/Type C K-Low Chemical			

For more information or samples, contact your distributor or visit: www.halyardhealth.com/cleanroom



¹ Not formulated with these commonly used vulcanizing chemicals: Sulfur, Thiurams, Thioxoles, Guanidines and Carbamates.
² Tested per ASTM D6319, EN 455-2
 *Registered Trademark or Trademark of O&M Halyard or its affiliates. ©2020, 2022. All rights reserved.
 H02757 H00078-12-01, 008271, COPY-04351 June 2022